



ThermoFisher
SCIENTIFIC



Raise the Bar

Environmental and Food Safety

Sofia, October 2018

Simonas Rudys
Technical Sales Manager
Emerging Markets EMEA
Chromatography and Mass Spectrometry

We Are The World Leader in Serving Science

Global Scale

- Approximately 70,000 employees globally
- More than \$20 billion in annual revenue
- Unparalleled commercial reach

Unmatched Depth

- Leading innovative technologies
- Deep applications expertise
- Premier laboratory productivity partner



ThermoFisher
S C I E N T I F I C

thermo
scientific

applied
biosystems

invitrogen

 fisher
scientific

unity
lab services

Our Mission: To enable our customers to make the world healthier, cleaner and safer

Our Business Segments

Analytical Instruments

Mass Spectrometry



Q-Exactive HF mass spectrometer

Chromatography



Vanquish UHPLC

Electron Microscopy



Titan S/TEM

Chemical Analysis



Gemini handheld analyzer

Life Science Solutions

Clinical Oncology



Next Gen Sequencing

Genetic Sciences



QuantStudio Dx R qPCR

Reproductive Health



Microarrays

Biosciences



Life Science Reagents

BioProduction



Cell Culture Reagents

Human Identification



GlobalFiler PCR Amplification Kit

Laboratory Products and Services

Lab Equipment



Sorvall WX+ Ultracentrifuge

Distribution and BioPharma Services



Enterprise and Instrument Services



Lab Consumables



E1 ClipTip Pipette system

Specialty Diagnostics

ImmunoDiagnostics



ImmunoCAP Allergy and EliA Autoimmunity Tests

Clinical Diagnostics



Oral-Eze Oral Fluid Collection System

Microbiology



Antimicrobial Susceptibility Testing Solutions

Pathology




Rotary Microtome

Transplant Diagnostics



NXType High Resolution Genotyping

Chemicals



A Laboratory Solutions Powerhouse



Software and Services

Comprehensive services
Multi-vendor support
Thermo Scientific™
Dionex™ Chromleon™
CDS, LIMS

Consumables

IC, LC, GC columns
Vials, septas, fittings
Solvents and standards



Instruments

GC, LC, IC
LC multiplexing
SQ, QQQ, Ion Trap,
Orbitrap, Hybrids

AA, ICP, ICP-MS, IRMS

**Expanded depth and scale offer
unequaled portfolio breadth**

Thermo Scientific LC Portfolio Overview

HPLC Systems

Routine Analysis

UHPLC Systems

Low-Flow

Flexibility, More Performance and Method Development

Research & HT

Vanquish Duo Workflows

Vanquish Horizon

- High pressure binary solvent mixing
- Industry-leading precision and accuracy
- Two thermostating modes
- Unmatched detection sensitivity
- Biocompatible
- Integration of multiple detection technologies

Vanquish Flex

- High pressure binary and low pressure quaternary solvent mixing options
- Two thermostating modes
- Biocompatible
- Integration of multiple detection technologies

UltiMate 3000 RS/BioRS

- Specialty workflow support
- Binary and Quaternary UHPLC systems

UltiMate 3000 RSLCnano

EASY-nLC 1200

- UHPLC systems for Nano/Cap/Micro range
- Ideal front-end for proteomics applications

UltiMate 3000 SD

- Workhorse for standard HPLC applications

UltiMate 3000 Basic Automated

- Highly economic & reliable



620 bar

800-1200 bar

Up to 1000 bar

1500 bar

Three Workflows. Two Flow Paths. One Integrated UHPLC Solution.




Dual LC

- A two channel LC system with standard footprint



Tandem LC/LC-MS

- System for efficient detector utilization



Inverse Gradient

- System for advanced quantification capabilities

Standard LC & LC-MS



Value

Standard instrumentation
Easy to use


Limitation

Idle time limits productivity
Low throughput
CAD quantification of unknowns

Target Market



Vanquish Duo for Dual LC




Value

Doubled sample throughput
Run complementary assays for better sample knowledge
Bench space savings

Target Market



Vanquish Duo for Tandem LC or LC-MS



Value

Increase LC-MS productivity by 50% with only 5% higher costs (LC-QQQ)

Target Market



Vanquish Duo for Inverse Gradient



Value

Standard-free quantitation of unknown compounds
Do not miss peaks

Target Market



Thermo Scientific Vanquish Duo System Benefits



- Faster from sample to result
- Higher sample throughput
- Efficient bench space utilization
- Improved return on investments
- Streamlined workflows by intelligent software
- Optimized quantification with Charged Aerosol Detector


Vanquish Duo

3 workflows requiring different instrument configurations

For all workflows



Workflow kits



Chromeleon Tools

e.g.wizards or method translation

Vanquish Duo for
Dual LC

Vanquish Duo for
Tandem LC/LC-MS

Vanquish Duo for
Inverse Gradient

Improved ROI by extended Productivity for each workflow

Thermo Scientific™ Vanquish™ Flex Dual Pump




- 2 ternary pumps in one housing
- 1000 bar (15,000 psi) version
- Delivers 2 completely independent gradients

Thermo Scientific™ Vanquish™ Dual Split Samplers




- 2 injection units in one autosampler
- 1000 and 1500 bar (15,000 or 22,000 psi) version
- Independent flow paths

Thermo Scientific™ Vanquish™ Duo System for Dual LC



Thermo Scientific Vanquish Duo for Dual LC – Instrument Configuration




What Can Be Accomplished with Thermo Scientific Vanquish Duo System for Dual LC?

Dual LC



Run two identical columns




Be finished in half the time

Finished in half the time

Dual LC



Run your complementary



Get more information out of every sample



Efficient characterization

Thermo Scientific Vanquish Duo for Tandem LC or LC-MS



Analytical Challenge – Thermo Scientific Vanquish Duo for Tandem LC or LC-MS

High investments in modern instrumentation



Standard LC / LC-MS




Time not acquiring meaningful data

Not fully utilizing instrumentation


Asked to maximize the instrument outcome

Thermo Scientific Vanquish Duo for Tandem LC or LC-MS – Instrument Configuration




Key Benefits – Thermo Scientific Vanquish Duo for Tandem LC or LC-MS


Advanced column wash without sacrificing throughput



Standard gradient profile




Tandem mode allows for time savings




Extra column washing steps without sacrificing throughput

Thermo Scientific Vanquish Duo for Inverse Gradient



Unique Detectors: Charged Aerosol Detection

- Detection of substances without a chromophore
- Consistent analyte response independent of chemical structure
- Dynamic range up to four orders of magnitude and highest sensitivity
- Relative quantification using virtually any standard




Comparison of Charged Aerosol Detection to UV and MS




Saving costs for additional detectors or assays

Thermo Scientific Vanquish Duo for Inverse Gradient – Instrument Configuration




Inverse Gradient Compensation for Uniform Response


Analytical Gradient



Analytical Gradient




Inverse Gradient




Eluent at detector

Eluent at detector

Conventional Gradient Elution



Inverse Gradient Compensation



Workflow Positioning



Improved Return on Investment

Saves bench space

Easy to use

Doubled sample throughput

Run complementary assays

Increased throughput for LC-MS


Reliably quantify unknown



Vanquish Duo
for Dual LC



Vanquish Duo for
Tandem LC or LC-MS



Vanquish Duo
for Inverse Gradient

✓ Double throughput or increase sample knowledge	✓ Increased detector utilization	✓ Reliable quantitation
✓	✓	✓
✓	✓	✓
✓	✗	✗
✓	✗	✗
✗	✓	✗
✗	✗	✓

New Sample Detection Portfolio

Integrated detector options

Advantages



Thermo Scientific™ Vanquish™ Horizon
Diode Array Detector

Outstanding signal to noise performance

NEW



Thermo Scientific™ Vanquish™ Flex
Diode Array Detector

Leading linearity, robustness and smooth MS
integration



Thermo Scientific™ Vanquish™ Flex
Variable Wavelength Detector

Supports a wide operating range



Thermo Scientific™ Vanquish™ Flex
Fluorescence Detector

High selectivity and sensitivity



Thermo Scientific™ Vanquish™ Horizon/Flex
Charged Aerosol Detector

Offers sensitive, universal detection with a nearly
uniform response

NEW



Thermo Scientific™ ISQ™ EC/EM Mass Spectrometer

Single quadrupole MS for robust and easy mass
detection

Three Workflows. Two Flow Paths. One Integrated UHPLC Solution.

Ultra
High
Productivity
Liquid
Chromatography




Dual LC

- A two channel LC system with standard footprint



Tandem LC or LC-MS

- System for efficient detector utilization



Inverse Gradient

- System for advanced quantification capabilities

Thermo Scientific Portfolio of Mass Spectrometry

- Metabolomics
- Proteomics
- Lipidomics
- Bioanalysis
- Clinical Research/
Precision Medicine



Non-targeted analysis

Orbitrap
HRAM



- Biomarker Discovery
- Proteomics
- Metabolomics
- Lipidomics
- Metabolism

Quantitative



- Food Safety
- Environmental
- Clinical Research/
Forensic Toxicology

Applied
Markets

For the
scientific
rigor you
demand

Research
Markets

Qualitative

- PTM Analysis
- Proteomics
- Forensic Toxicology




MS, MSⁿ


Targeted analysis

2018 LSMS Hardware Portfolio

Triple Quadrupoles



TSQ Altis



TSQ Quantis



TSQ Fortis



Q Exactives



Q Exactive HF-X



Q Exactive HF




Q Exactive (Plus)




Q Exactive Focus

Hybrids



**Options:
UVPD, 1M, APD**



**Orbitrap Fusion
Lumos**



Orbitrap Fusion



Orbitrap Elite

Confident quantitation for any compound, any matrix, any user

PERFORMANCE

Environmental and Food Safety
Clinical Research
Pharma QA/QC



TSQ Fortis

- Mass Range m/z 5 – 3000
- Max Resolution **0.4 FWHM**
- Max 30,000 transitions per run
- Polarity Switching < 20 msec
- Dynamic interscan time
- 600 SRM/sec
- TNG software
- Chromeleon support
- **50,000:1 S/N**

Food Safety
Pharma
Clinical Research
Forensic Toxicology



TSQ Quantis

- Mass Range m/z 5 – 3000
- Max Resolution **0.4 FWHM**
- Max 30,000 transitions per run
- Polarity Switching < 20 msec
- Dynamic interscan time
- 600 SRM/sec
- TNG software
- Chromeleon support
- **150,000:1 S/N**

Pharma/Biopharma
Environmental and Food Safety
Omics



TSQ Altis

- Mass Range m/z 5 – 2000
- Max Resolution **0.2 FWHM**
- Max 30,000 transitions per run
- Polarity Switching < 20 msec
- Dynamic interscan time
- 600 SRM/sec
- TNG software
- Chromeleon support
- **500,000:1 S/N**


VALUE

TSQ Fortis: Affordable Productivity, For Everyone



Active Ion Management Plus (AIM+) - The next step in precision design delivers the ultimate in ion management, inception to detection, from the OptaMax™ ion source housing to the enhanced electron multiplier. Incorporates segmented quadrupoles with hyperbolic surfaces and enhanced RF Electronics to further optimize ion management precision, reliability, speed, and reproducibility.

Enhanced dual-mode electron multiplier detector ensures excellent linearity and dynamic range



Segmented Quadrupoles with hyperbolic surfaces for enhanced performance with both SRM and H-SRM (0.4 FWHM)



Active collision cell with axial DC field facilitates more SRMs/sec

Matrix Separator Ion Guide (MSIG) Ensures Robustness while efficiently transmitting the ion beam



Ion beam guide with neutral blocker Reduces chemical background



OptaMax™ NG
APCI ready

Matrix Separator Ion Guide (MSIG) – Unique to TSQ Fortis!



Benefits: Increased Robustness while efficiently transmitting the ion beam and ensuring the required sensitivity

- A mass dependent potential is applied to focus the ions from the ion guide region toward the opening of the matrix separator.
- The off-axis design of the matrix separator enables the separation of neutrals from the ion beam, leading to increased robustness and therefore increased uptime.

Segmented Quadrupoles

Benefits: Increased Sensitivity, Flat tuning for consistent and robust performance

- The use of RF only pre-filters (segments) between the entrance lens and the quadrupole minimizes the effects of fringe fields, leading to improved transmission (and therefore sensitivity) at unit and higher resolution.
- With the RF only pre-filter, the tuning of several lenses is flat across mass range allowing the voltage to be set and not tuned. This helps reducing the complexity of the tune and making the systems more consistent.



*New segmented quadrupoles
with hyperbolic surfaces*

Round Rods vs Hyperbolic rods





Competition



Thermo Scientific

Effect of Peak Width On Transmission




TSQ Altis – 0.2 Da FWHM

TSQ Quantis – 0.4 Da FWHM

TSQ Fortis – 0.4 Da FWHM


Hyperquads, Precursor Isolation and H-SRM

- Enhanced isolation and stability of ion(s)
- H-SRM – Improved isolation of precursor ion(s)
- More Selectivity – isolate only ion(s) of interest
- **Benefit** - Improved signal to noise
- **Benefit** More Confidence – less false positives
- **Benefit** Higher Quality data



Improved isolation of precursor without significant loss of signal

Clenbuterol Urine Extract



Unit Resolution 0.7 Da FWHM

H-SRM 0.2 Da FWHM




Figure 1. SRM vs. H-SRM for Chloramphenicol (321→152) of a 0.03 µg/kg fortified pork kidney sample.

What Makes the new Triple Quads Robust?



Robustness

The ability to perform at the expected level (LOD, LOQ, MRL) under adverse conditions (complex samples, limited sample preparation) for the desired period of time without maintenance.

OptaMax NG



Source housing

Reliable and consistent performance with improved usability!


ITT/Sweep Cone



Ion Transfer Tube and Sweep Cone

No need to break vacuum for basic maintenance!

MSIG




Matrix Separator Ion Guide

TSQ Fortis only!

Eliminates neutrals while ensuring efficient transmission and sensitivity


Ion Beam Guide



Ion beam guide with neutral blocker

Blocks neutrals - keeps the ion path clean and reduces chemical background!

What Makes the new Triple Quads Consistent?



Consistency

The ability to report equivalent results every time independently of the system or the user (day-to-day and system-to-system)


OptaMax NG



Re-designed Source Housing

Reliable and consistent performance with improved usability!


Optics



Improved Optics Alignment

More consistent performance from system-to-system!

Segmented Quads



Segmented Quads with hyperbolic surfaces

Flat tuning across the mass range!

TSQ Fortis: Demonstration of Robustness – Clinical Research

Immunosuppressants





+2000 injections over 6 days
(CV% ≤ 6%)

Red lines represent ± 20% of
calculated amounts (ng/mL)

Yellow line represent the period in
which the ion transfer tube was
cleaned (user basic maintenance
~5 min operation) to demonstrate
consistent performance before
and after user maintenance.




Technical Note 65206




TSQ Quantis: Demonstration of Robustness – Food Safety



Atrazine QC monitored in leek for more than 400 injections with 4.5% RSD. Red lines represent $\pm 20\%$ response at $10 \mu\text{g}/\text{Kg}$. Yellow lines show the time the system was placed in standby mode for 12h to demonstrate consistent performance after standby period.




Application Note 64971



Thermo Scientific
Pesticide Explorer Collection


Start-To-Finish Workflows Pesticide Analysis

Thermo Scientific Pesticide Explorer Collection for Quantification



Thermo SCIENTIFIC





Multi-class Veterinary Drug Screening and Quantitation with a Comprehensive Workflow

Balancing the Complexity of Multi-Class Veterinary Drug Method Development


Several single-class methods



One multi-class method




Stability of multi-component mixtures



Generic sample prep with good recoveries + RSDs




Wide range of chemical classes+ MRLs



Inert LC system and column for wide pKa range and good peak shape

Components in the Workflow


QuEChERS
Sample Prep
& Columns




UHPLC



Sensitive
QQQ



Complete
Methods with SRM,
compound database



**QuEChERS
Sample
Preparation**

(Consumables +
procedure)

**UHPLC Chromatographic
Separation**

(Inert LC pathway + Column)

QQQ Detection

(Robust system with
sensitivity to meet global
MRLs over diverse cmpd
classes)






**Data Processing
&
Reporting**

Technical Highlights of a Comprehensive Workflow

Feature	Detail	Comment
# of Compounds	170+ Analytes with acquisition and master processing method	Compound data base with optimized SRMs for easy set-up
Matrices	Muscle meat (cattle), milk, and salmon (fillet)	Demonstrates broad applicability of the multi-class method approach
QC and Analytical Standards	QC Check-20 cmpds Standards- 170 compounds	Ability to analyze 'System Performance Check'
Sample Preparation	Generic QuEChERs extraction w/simple clean-up	Detailed procedure required for applicable matrices
Column	Thermo Scientific™ Accucore™ VDX 100mm x 2.1 x 2.6 um	Ensure performance across wide polarity range
UHPLC-MS system	Thermo Scientific™ TSQ Altis™ Triple Quadrupole MS Thermo Scientific™ Vanquish™ Flex Binary pump	Includes detailed User guide for proper installation and system check-out
Software	Thermo Scientific™ TraceFinder™ Software 4.1 SP5	Unify all aspects of data handling from sample acquisition to report


Chemical Classes within the Method

Compound Classes



Sample Preparation and LC Conditions

- **QuEChERS based approach**
 - EDTA/NH₄ oxalate solution and acetonitrile
 - Sample homogenised until fully dispersed
 - Sodium sulphate added before centrifugation
 - Dispersive SPE (CEC-C₁₈) clean-up
 - Add 1 mL H₂O to 3mL extract, filter, inject
- **LC conditions**
 - Thermo Scientific™ Accucore™ VDX Column
2.1 x 100 x 2.6 μm
 - MP A: 0.05% formic acid
 - MP B: 0.05% formic acid, 5% H₂O in 1:1 MeOH:MeCN
 - 2 μL injection
- **Acquire Data on TSQ Altis**
 - *Use pos/neg switching*
 - Comprehensive CDB with all optimized SRMs






Column-Accucore VDX- Robust and Selective for VetDrugs




Features

- Solid core particle- high resolution separations
- Column chemistry selectivity similar to C18 columns
- Optimized for MS detection
- Low column bleed
- Optimized for low tailing
- **Robust against matrix extracts**
- Particle size: 2.6 μm



Improved Sensitivity and Specificity with H-SRM (0.2 Da FWHM)




SRM, 0.7 FWHM



H-SRM, 0.2 FWHM

Advantage of H-SRM for Fluazuron in salmon fillet extract: Noise is significantly reduced allowing improved signal-to-noise (selectivity) against the matrix with high ion transmission




Calibration curve for Flunixin in salmon fillet matrix using H-SRM.


Compound	MDL (0.7 SRM)	MDL (0.2 HSRM)
Derquantel	0.42	0.04
Fenbendazole Sulfone	0.33	0.05
Flunixin	0.38	0.04
Hydroxyipronidazole	0.62	0.08
Robenidine	0.33	0.05
Teflubenzuron	4.47	0.47

Advantage of H-SRM for improved method detection limits (MDLs) for select compounds in a salmon fillet extract

TIC and SRMs




TIC of all SRMs @ 1xSTC in TraceFinder software




The distribution of 546 SRMs are easily visualized in the instrument software tune page

Quantitative Results- 0.2 to 5 x STC-Bovine



Sulfamethazine in bovine extract at 1× STC, with screening range from 10-250 ng/g.


Quantitative Results- 0.2 to 5 x STC - Salmon Fillet




Ethyl Violet in salmon extract @ 1×STC, with a screening range from 0.2–5 ng/g.

Compound Class - Average Calculated MDL (ng/g)

Average Method Detection Limit (MDL) by Class



Steps for Evaluating Method Performance



Validate precision on 4 separate LC/MS/MS Systems to ensure robustness and method transferability

Calculate Absolute % Recovery Based upon a 'post spike' at 3 x STC

Calculate MDL @ or below a cut-off (Lowest xSTC factor at or below 15% RSD)

Calculate %RSDs at each level to check precision

Prepare Matrix Extracted Spike (MES) Replicates at 0.2 - 5 x STC (Establish 'Calibration Line' for screening)

Establish Screen Target Concentration (STC) Level - Levels typically 1/4 to 1/3 the MRL

Retention Time Stability Across 4 LC/MS/MS Systems


- Validation included 4 different LC/MS/MS systems: 3 Thermo San Jose Labs and 1 at Iowa State Lab
- Representative compounds throughout the gradient RT range show excellent stability in matrix (Salmon Fillet shown), with excellent agreement between different systems (different columns, 3 LCs, and MS systems).

Compound	Exp RT	RT System 1	RT System 2	RT System 3	RT System 4	AVE	SDEV	%RSD
Cyromazine	0.83	0.82	0.82	0.85	0.82	0.83	0.0150	1.81
Dimetridazole	4.16	4.07	4.11	4.15	4.12	4.11	0.0330	0.80
Sulfamethoxazole	5.84	5.80	5.81	5.86	5.79	5.82	0.0311	0.53
Oxyphenbutazone	8.52	8.49	8.49	8.53	8.49	8.50	0.0200	0.24
Triclabendazole	10.21	10.19	10.22	10.24	10.21	10.22	0.0208	0.20
Closantel	11.88	11.91	11.91	11.92	11.93	11.92	0.0096	0.08
Monensin	12.21	12.23	12.25	12.27	12.26	12.25	0.0171	0.14

Example Data : Salmon Matrix- Precision at STC Levels Across 4 LC/MS/MS Systems

	Lowest Global								
Compound	MRL	STC	System 1	System 2	System 3	System 4	Average	SDEV	%RSD
Amoxicillin	50	25	24.4	24.8	25.0	27.3	25.4	1.31	5.2
Chlortetracycline	200	100	87.3	97.1	86.4	104.0	93.7	8.39	9.0
Danofloxacin	100	50	50.9	49.8	47.3	52.3	50.1	2.09	4.2
Dicloxacillin	300	150	144.5	148.2	147.7	152.2	148.2	3.15	2.1
Difloxacin	300	150	150.0	149.2	145.6	159.3	151.0	5.80	3.8
Doxycycline	10	10	9.2	9.1	9.2	10.1	9.4	0.47	5.0
Enrofloxacin	100	50	50.2	49.4	55.0	52.3	51.7	2.51	4.9
Flumequine	500	200	200.8	201.9	201.5	199.3	200.9	1.17	0.6
Oxolinic Acid	50	25	28.2	29.3	27.5	25.8	27.7	1.48	5.3
Oxytetracycline	200	100	102.2	96.4	94.6	101.1	98.6	3.65	3.7
Penicillin G	50	25	23.6	24.0	21.6	25.4	23.6	1.57	6.7
Sarafloxacin	10	3	3.0	3.0	2.9	3.1	3.0	0.06	2.0
Sulfadoxine	SUM 100	10	9.1	8.7	9.4	10.1	9.3	0.59	6.3
Tetracycline	200	100	98.8	98.3	91.7	100.4	97.3	3.81	3.9
Thiamphenicol	50	25	23.6	24.0	22.2	26.5	24.1	1.80	7.5
Tilmicosin	50	25	27.5	28.0	25.4	23.2	26.0	2.21	8.5
Trimethoprim	50	25	21.8	18.7	24.0	22.0	21.6	2.19	10.1
Tylosin	100	50	46.8	47.5	45.9	47.7	47.0	0.80	1.7

Chromatography with mass spectrometric detection



Current golden standard:


Targeted quantitative measurement by LC-MS/MS and GC-MS/MS

Multi-residue methods, ~150–250 analytes/method

Emerging:

Non-targeted measurement by LC and GC + full scan MS
for better coverage of pesticide scope and easier measurement
IC for easier separation of polar substances

Sample preparation/analyte detection strategy



Former Pesticide Multi-Residue Method Setup




■ Extraction

Acetonitrile, Ethyl acetate, Methanol...



■ Clean-up

GPC, SPE, LLE, LC



■ Determination

GC, LC, GC-MS, LC-MS, GC-MS/MS, LC-MS/MS...

Mostly replaced by **QuEChERS** today



Thermo Scientific™ QuEChERS™ method

Improving QuEChERS extraction tips & tricks:

- **Dry food (cereals/dried food, < 25 % water content):**
 - addition of water to enable adequate partitioning and reducing interaction of pesticides with matrix
- **Food containing fat/wax (avocado/oil):**
 - after extraction step add a freezing out step and transfer supernatant to clean-up tube
 - more clean-up might be needed of raw extract (PSA+C18)
- **Food containing complex matrix (tea/spices)**
 - additional clean-up with GCB might be necessary (potential loss of planar structure pesticides like thiabendazole)
- **Acidic food (citrus):**
 - adjust pH (5-5.5) to increase recovery (e.g. citrate buffering salts in QuEChERS extraction tube) and reduce coextraction of matrix interferences (Note: acid labile compounds require higher pH 8)



Improving QuEChERS Clean-up Tips & Tricks:


- **QuEChERS clean-up tube additives:**

Product Selection

Matrix Type	Examples	Sorbent Requirements
General Matrices	Apples Cucumber Melon	Magnesium Sulfate, PSA
Fatty Matrices	Milk Cereals Fish	Magnesium Sulfate, PSA, C18
Pigmented Matrices	Lettuce Carrot Wine	Magnesium Sulfate, PSA, C18, GCB
High Pigmented Matrices	Spinach Red Peppers	Magnesium Sulfate, PSA, C18, GCB

[Click here for more information](#)

What is the future?



Various Screening Approaches | Common Definitions

	Targeted Screening	Untargeted Screening	Unknown Screening
SEARCH FOR:	what you know	what you suspect	all the rest
HOW YOU SEARCH:	Compound search after targeted acquisition; uses specific compound list	Compound search follows acquisition based on most intense signals NOT on a list of compounds	Search of compounds with a totally untargeted acquisition
ADVISORY:	Detection parameters need to be numerous & selective to avoid false positives	Identification based on compound databases and/or spectral libraries	Putative identification typically followed by structural elucidation

Exactive Series Product Portfolio

Increased Performance

- Orbitrap analyzer
- Mass Range m/z 50 - 6000
- Mass Range EMR: 300 - 20,000
- Mass Accuracy: <1ppm
- Mass Resolution >140,000
- Scan Speed up to 12 Hz
- Octapole transfer for MS and AIF



Exactive Plus (EMR)

- Orbitrap analyzer
- Mass Range m/z 50 - 3000
- Mass Accuracy <1ppm
- Max. Mass Resolution >70,000
- Scan speed up to 12Hz MS and MS/MS
- Polarity switching
- Hyperbolic quadrupole: MS, ddHCD, AIF, SIM, PRM, vDIA
- Multiplexing SIM only
- Refined workflows for applied markets



Q Exactive Focus

- Orbitrap analyzer
- Mass Range m/z 50 - 6000
- Mass Accuracy <1ppm
- Max. Mass Resolution >140,000
- Scan speed up to 12Hz MS and MS/MS
- Polarity switching
- Hyperbolic quadrupole: MS, ddHCD, AIF, SIM, PRM, DIA
- Spectral Multiplexing



Q Exactive



Q Exactive Plus

- Orbitrap analyzer
- Mass Range m/z 50 - 6000
- Mass Accuracy <1ppm
- Max. Mass Resolution >140,000
- Scan speed up to 12Hz MS and MS/MS
- Advanced Quadrupole Technology (AQT)
- Advanced Active Beam Guide (AABG)
- Polarity switching
- Segmented quadrupole: MS, ddHCD, AIF, SIM, PRM, DIA
- Spectral Multiplexing
- PRM dyn-RT correction
- Opt. Intact Protein Mode
- Opt. Enh Res. Mode (280k)
- Optional BioPharma Mode mass range up to 8000



Q Exactive HF


- Ultra High Field Orbitrap analyzer
- Mass Range m/z 50 - 6000
- Mass Accuracy <1ppm
- Max. Mass Resolution >240,000
- Scan speed up to 18Hz MS and MS/MS
- Advanced Quadrupole Technology (AQT)
- Advanced Active Beam Guide (AABG)
- Polarity switching
- Segmented quadrupole: MS, ddHCD, AIF, SIM, PRM, DIA
- Spectral Multiplexing
- PRM dynamic RT correction
- Opt. Intact Protein Mode
- Optional BioPharma Mode mass range up to 8000



Q Exactive HF-X

- Ultra High Field Orbitrap analyzer
- Mass Range m/z 50 - 6000
- Mass Accuracy <1ppm
- Max. Mass Resolution >240,000
- Scan speed up to 40Hz MS and MS/MS
- High capacity transfer tube
- Electrodynamic ion funnel
- Advanced Quadrupole Technology (AQT)
- Advanced Active Beam Guide (AABG)
- Polarity switching
- Segmented quadrupole: MS, ddHCD, AIF, SIM, PRM, DIA
- Spectral Multiplexing
- PRM dynamic RT correction
- Advanced ddHCD algorithm
- Optional BioPharma Mode mass range up to 8000 with Intact Protein settings

Schematic of the Q Exactive Focus



Benefits of Using Orbitrap Technology



Resolving Power

Separate target compounds from interference



Mass stability and mass accuracy

Calibration lasts days to weeks because of unrivaled scan:scan mass accuracy & precision



Fast polarity switching

Quick positive-negative switch maintaining accurate mass/charge determination



No trade-off between resolution and sensitivity

Both resolution and sensitivity are retained



Minimum scan to scan variance

Excellent signal/noise, no need for averaging






Experiment flexibility

Easy to use, a plug and play device


All your screening and quantitation applications in One platform

Selectivity increases with higher mass accuracy



*Zachariasova M et al, Anal. Chim. Acta

Does resolution matter?



RAFA 2013: Łukasz Rajski, María del Mar Gómez Ramos, Amadeo R. Fernández-Alba; EURL for Pesticide Residues in Fruits and Vegetables. Pesticide Residue Research Group. University of Almería, Spain. e-mail: amadeo@ual.es

3 ways of Quantitation/Screening for Routine Work

Full MS or targeted SIM/ddMS2

- Post-acquisition - extracted ion chromatograms of parent ions of interest
- Relies on high resolution for selectivity
- Useful for less complex background
- No method development/preparation needed

Full MS/ All Ion Fragmentation – vDIA*

- Post-acquisition - extracted ion chromatograms of parent ions of interest
- Scheduled target (inclusion) list (Rt, m/z)
- Minimum method development (e.g., predefine parent ions, tr)
- Also for screening purposes

PRM (Parallel Reaction Monitoring)

- Post-acquisition – extracted ion chromatograms of parent -> fragment transitions acquired
- Scheduled target list (Rt, m/z , collision energy)
- Most sensitive and selective even in highly complex matrices

Experiments

General

- Full MS
- SIM
- PRM
- Full MS - AIF
- Full MS - vDIA

Experiments


General

- Full MS
- SIM
- PRM
- Full MS - AIF
- Full MS - vDIA


Experiments

General


- Full MS
- SIM
- PRM
- Full MS - AIF
- Full MS - vDIA



Q Exactive Focus Scan Methods



Orbitrap MS: Unmatched Resolution



Benefits of high resolution for small molecule analysis

- Differentiate similar masses in complex matrix
- Isobaric species
- Fine isotopic pattern

High resolution and mass accuracy increase confidence of unknown identification

Analysis of Polar Pesticides by IC-MS/MS

IC vs LC arguments

	Ion Chromatography		Liquid Chromatography	
Matrix	Food	Water	Food	Water
Sample preparation	Extraction (QuPPE)	None (filtration)	Extraction (QuPPE) (optional FMOCC* derivatization)	FMOCC* Derivatization (SPE cleanup) (automation possible)
Separation quality	Excellent	Excellent	Poor (Hypercarb)	Very good (Reversed Phase)
LOQs	1-10 ppb	5 -10 ppt	1 – 10 ppb	5 – 10 ppt
MS/MS	**TSQ Endura/Quantis	**TSQ Quantiva/Altis	**TSQ Quantiv/Altis (Endura not suitable)	**TSQ Quantiva/Altis
Matrix tolerance	Very good	Very good	Limited	Good
Time consumption	Low	Low	Low - medium	Medium – High (automation possible)



*FMOCC derivatization applicable only for Gly, AMPA, Glu!

**Thermo Scientific™ TSQ™ Endura/Quantis, Thermo Scientific™ TSQ™ Quantiva/Altis™




TSQ Fortis (& Ion Chromatography): Environmental Analysis

Haloacetic Acids by IC-MS/MS

Compound	MDL (µg/L)
MCAA	0.03
MBAA	0.03
DCAA	0.02
DBAA	0.02
BCAA	0.03
TCAA	0.06
BDCAA	0.05
DBCAA	0.15
TBAA	0.15
Dalapon	0.03
Bromate	0.02



Perchlorate by IC-MS/MS



Application Note 65196

Application Note 65201


IC-MS/MS solutions for environmental analysis Quantitation of **Glyphosate** at 5 ng/L, 3% RSD

Experimental Details

ICS 5000+
Flow rate: 0.3 mL/min
Eluent Source: Eluent Generator
Mobile Phase: KOH

MS: TSQ Altis

Software: TraceFinder Software 4.1




Polar pesticides in fruits and vegetables by IC-MS/MS Instrumental Configuration

- IC-System: Dionex Integrion HPIC
- Eluent Source: Dionex EGC 500 KOH
- Suppressor: ASRS 300
- Autosampler: AS – AP Dionex Autosampler
- Mass Spectrometer: TSQ Quantiva

- Column: AS24 (2 x 250 mm)
- Guard Column: AG24 (2 x 50 mm)
- Eluent: KOH
- Column Temperature: 30 °C
- Flow rate: 0.3 ml/min
- Make-up flow: 0.1 ml/min
- Make-up solvent: Isopropanol (0.1 mL/min)
- Run time: 22 min
- Injection volume: 10 µL
- Gradient KOH:

No	Time	Concentration [mM]	Curve
1	{Initial Tim	Instrument Setup	
2	{Initial Tim	22.00	5
3	0.000	Run	
4	0.000	22.00	5
5	4.100	22.00	5
6	7.000	25.00	5
7	7.100	40.00	5
8	9.500	60.00	5
9	12.000	80.00	5
10	14.500	80.00	5
11	15.000	100.00	5
12	17.000	100.00	5
13	17.100	22.00	5
14	19.900	22.00	5
15	20.000	22.00	5
16	22.000	22.00	5
17	New Row		
18	22.000	Stop Run	



Compound	LOD (µg/kg)	LOQ (µg/kg)
Phosphonic acid	0.50	1.0
Fosetyl	0.03	0.1
AMPA	0.50	1.0
Glufosinate	0.50	1.0
N-acetyl-glufosinate	0.25	0.75
Ethephon-hydroxy	0.25	0.75
N-acetyl-AMPA	0.10	0.50
MPPA	0.10	0.50
Ethephon	0.25	0.75
Glyphosate	0.03	0.10



Compound	LOD (µg/kg)	LOQ (µg/kg)
Phosphonic acid	0.50	1.00
Fosetyl	0.03	0.10
AMPA	0.50	1.00
Glufosinate	0.25	0.75
N-acetyl-glufosinate	0.25	0.75
Ethephon-hydroxy	0.03	0.10
N-acetyl-AMPA	0.03	0.10
MPPA	0.03	0.10
Ethephon	0.25	0.75
Glyphosate	0.03	0.10



Fast Analysis of Polar Pesticides in Water by IC-MS/MS

IC: Dionex Integrion™
MS: TSQ Quantiva™



Column: AS24 (2 x 250 mm)
Guard Column: AG24 (2 x 50 mm)
Eluent: KOH
Injection volume: 100 µL
Column Temperature: 21 °C
Flow rate: 0.3 ml/min
Make-up flow: 0.1 ml/min
Make-up solvent: CH₃OH

Time (min)	Concentration of KOH in eluent (mM)
0	25
0.2	25
11	80
11.1	100
12.5	100
12.6	25
17.0	25

Component	Matrix	LOD [ppt]	LOQ [ppt]
Fosetyl	Drinking	2.5	5
	Bottled	1	2.5
	Surface	2.5	5
Glufosinate	Drinking	5	10
	Bottled	5	10
	Surface	5	10
AMPA	Drinking	5	10
	Bottled	5	10
	Surface	5	10
Clopyralid	Drinking	10	50
	Bottled	5	10
	Surface	5	10
Glyphosate	Drinking	10	50
	Bottled	10	50
	Surface	10	50

Thermo Scientific Products in Food Testing

- Sample preparation
- Chromatography
- Mass Spectrometry
- Trace Elemental
- IRMS
- Consumables
- Services
- Informatics

