

**ThermoFisher**  
SCIENTIFIC


## Dioxin and PCBs with Triple Quadrupole and HRMS Technology- Setting the New Standard

Michal Godula  
Special Solutions Center  
Thermo Fisher Scientific

The world leader in serving science


# Persistence & Accumulation in the Food Chain

**POPs are chemically stable environmental contaminants**



## UNEP definition of Persistent Organic Pollutants (POPs)

“...chemical substances that **PERSIST** in the environment, **BIOACCUMULATE** through the food web, and pose a **RISK** of causing adverse effects to **HUMAN HEALTH** and the environment. With the evidence of long-range transport of these substances to regions where they have never been used or produced and the consequent threats they pose to the environment of the whole globe, the international community has now, at several occasions called for urgent global actions to **REDUCE AND ELIMINATE RELEASES** of these chemicals.”



## Dioxins at Trace Levels are Highly Toxic

- Tolerable Daily Intake    1 pg/kg/d TEQ
- Amounts typically detected:
  - Food                            PPT
  - Water                            PPQ
  - Human Serum                PPQ
  - Air                                fg/m<sup>3</sup>


# Toxicity



- Viktor Yushchenko, President of Ukraine (2005-2010)
- Poisoned in 2004

## Basic Concepts of Toxic Equivalency Factor (TEF)

- Dioxin-like toxicity is **cumulative**
  - So we can compare toxicology of various dioxin-like compounds and simplify risk assessment and regulatory control.
- All other compounds with dioxin-like toxicity are compared to 2378 TCDD and given a relative TEF



- Toxic Equivalency Factor (TEF) of 2378-TCDD = 1

# Ever Changing TEF Concept

- TEF factors are continuously revised
- TEFs are species dependent
  - Human
  - Fish
  - Avian
  - Mammalian
- Generally we are interested in human exposure

Congener	I-TEF	WHO1998-TEF	WHO2005-TEF
<b>Polychlorinated dibenzo-p-dioxins</b>			
2378-CI4DD	1	1	1
12378-CI5DD	0.5	1	1
123478-CI6DD	0.1	0.1	0.1
123678-CI6DD	0.1	0.1	0.1
123789-CI6DD	0.1	0.1	0.1
1234678-CI7DD	0.01	0.01	0.01
CI8DD	0.001	0.0001	0.0003
<b>Polychlorinated dibenzofurans</b>			
2378-CI4DF	0.1	0.1	0.1
12378-CI5DF	0.05	0.05	0.03
23478-CI5DF	0.5	0.5	0.3
123478-CI6DF	0.1	0.1	0.1
123678-CI6DF	0.1	0.1	0.1
123789-CI6DF	0.1	0.1	0.1
234678-CI6DF	0.1	0.1	0.1
1234678-CI7DF	0.01	0.01	0.01
1234789-CI7DF	0.01	0.01	0.01
CI8DF	0.001	0.0001	0.0003
<b>Non-ortho polychlorinated biphenyls</b>			
PCB-77	-	0.0001	0.0001
PCB-81	-	0.0001	0.0003
PCB-126	-	0.1	0.1
PCB-169	-	0.01	0.03
<b>Mono-ortho polychlorinated biphenyls</b>			
PCB-105	-	0.0001	0.0003
PCB 114	-	0.0005	0.0003
PCB-118	-	0.0001	0.0003
PCB-123	-	0.0001	0.0003
PCB-156	-	0.0005	0.0003
PCB-157	-	0.0005	0.0003
PCB-167	-	0.00001	0.0003
PCB-189	-	0.0001	0.0003

- EU guidelines for food analysis use WHO2005 TEF factors

# Persistent Organic Pollutants

- Incineration products

- Polychlorinated dibenzodioxins –PCDDs (75)
- Polychlorinated dibenzofurans - PCDFs (135)

- Environmental contaminants

- Polychlorinated biphenyls – PCBs (209)
- Polybrominated diphenyl ethers – PBDEs (209)
- Chlorinated paraffins
- Organochlorine pesticides
- Polyfluorinated ether sulphones




Large numbers of congeners but of differing toxicity

High specificity is required to focus on toxicologically significant POPs




# Structural Similarities Requires High Specificity


## Polychlorodibenzo-p-dioxin





## Polychlorodibenzofuran



## Chlorinated Paraffins



## Polychlorinated biphenyls



# EU Regulations for POPs in Food & Feed

- Regulation (EC) No 1881/2006

For example:

- Sum of dioxins (WHO-PCDD/PCDF-TEQ)
  - Limits range from 0.75 pg/g fat for vegetable oils to 4.0 pg/g for fish
- Sum of dioxins & dioxin-like PCBs
  - Limits range from 1.5 pg/g fat for vegetable oils to 8.0 pg/g for fish
- Separate limits apply to animal feed

- Regulation (EC) No 396/2005

For example:-

- DDT (sum of *pp'*-DDT, *op'*-DDT, *pp'*-DDT and *pp'*-TDE (DDD) expressed as DDT
  - Limits range from 0.04 mg/kg for milk 0.05 mg/kg for fruit to 1.0 mg/kg for meat and offal

# “Dioxin Incidents” – Economic Impacts

## Dioxin and PCB incidents are expensive.

- Egg production enterprise with \$5M hens would cost in excess of \$30M.
- Broiler enterprise producing 3M broilers per week would exceed \$85M.
- Times Beach, Missouri, road dust covered by oil spread. The cleanup cost = \$110M million.
- The total costs of the Belgian food crisis are estimated up to \$1000M.



## Repeated “Dioxin” cases from only recent years:

- 1999 Belgian PCB/dioxin in eggs, poultry
- 1999 Clay and zeolithes for feed
- 2000 Choline chloride
- 2004 Potato pulp
- 2005 Hydrochloric acid
- 2007 Indian Guar Gum thickener
- 2008 Irish Pork Meat
- 2008 Italian Mozzarella Cheese ...
  - ... still on today around Naples.
- 2011 German eggs, poultry contaminated from feed



# Types of Matrices Requiring POP Analysis

- Agricultural products & foodstuffs
  - Animal & fish products, oils & fats
- Animal feedingstuffs
  - Compound feeds and feed components
- Environmental samples
  - Incinerator fly ash, soil, grass and marine biota
- Biological samples
  - Human and animal blood, human milk & adipose tissue




# How To Analyze Dioxins - Complete POPs Solutions

Step	Requirement	Solution
Laboratory PPE	Fume hoods, lab coats, gloves, eye protection etc.	Fisher Scientific
Laboratory Reagents	Alumina, Silica, Solvents, Acids	Fisher Scientific
Laboratory Apparatus	Freeze driers, ovens, balances, pipettes, vials, columns	Fisher Scientific
Sample Preparation	Extraction	Thermo Scientific™ Dionex™ ASE350 Accelerated Solvent Extraction
	Extract Fractionation / Purification	CAPE Technologies Sample Preparation Kits*
	Extract Concentration	Thermo Scientific™ Dionex™ Rocket™ Evaporator
Sample Analysis	Chromatography	Thermo Scientific™ Trace 1310 GC, TriPlus RSH & Trace Gold GC/MS columns and consumables
Sample Analysis	GC-HRMS or GC-MS/MS	Thermo Scientific™ DFS™ GC-HRMS or Thermo Scientific™ TSQ™ 8000 Evo GC-MS/MS
Data Handling & Reporting	Dioxin Data processing, QA/QC checking & reporting	Thermo Scientific™ TargetQua™ n 3 Thermo Scientific™ Xcalibur™
	LIMS	Thermo Scientific™ Nautilus™ LIMS

- Thermo Fisher Scientific supports the entire dioxin analysis workflow

# Typical Dioxin Workflow (EPA1613)



## Background Info: Confirmation vs Screening

- **Confirmatory Analysis**

- Unequivocal identification and quantification
- Results can be used in prosecution

- **Screening Analysis**

- Indication of levels only
- Must be followed with confirmatory analysis



# What's New In EU Regulation For Dioxins?

## Past

[Commission Regulation \(EU\) No 252/2012 of 21 March 2012](#)

- Specifies use of **GC-HRMS** for confirmatory dioxin analysis
- **GC-MS/MS** was allowed a screening technique.



## Present

[Commission Regulation \(EU\) No 589/2014 of 2 June 2014](#)

- Specifies use of **GC-HRMS** or **GC-MS/MS** for confirmatory dioxin analysis
- **GC-MS/MS** is “an appropriate confirmatory method for checking compliance with the maximum”, only.
- **GC-HRMS** remains the recommended technique for “determination of low background levels in food monitoring, following of time trends, exposure assessment of the population”.




# Definitions of Screening and Confirmation

## Screening

- High throughput – fast
- Low cost
- Low false negative rate (< 1%)
- False positives acceptable (<5%)
- Relaxed precision (<30 % RSD)
- Sensitivity 25% below of ML

## Confirmation

- Unequivocal identification
- High precision (<15 % RSD)
- Higher sensitivity to monitor lower background and action levels.
- High resolution necessary



Source: EU Directives 96/23/EC, 1883/2006


# Options for screening POPs in Food & Biological Samples

- **DR-CALUX**


- Receptor-based assay
- Gives indication of total toxicity of extracts
- Requires adequate sample concentration & cleanup
- Lacks flexibility – only PCDDs/PCDFs and PCBs determined

- **GC-MS/MS (also for confirmation)**

- Instrumental
- Can measure individual congeners & TEQs
- Requires adequate sample concentration & cleanup
- High flexibility – can be used for all POPs as well as other residues & contaminants




# Extraction & Cleanup for Dioxins & PCBs – HR-GC-MS



Time to response: 3-5 days

# Extraction & Cleanup for Dioxins & PCBs – GC-MS/MS



Provided by G. Brambilla, ISS Rome, Italy

# GC-MS/MS vs. GC-HRMS



**TSQ 8000 Evo**

**DFS**

GC-MS/MS

GC-HRMS

Triple Quadrupole MS

Magnetic Sector HRMS

High performance, easy to use MS/MS for non-experts

High-Resolution Full Scan and SIM

Target Analysis in Complex Matrix

Dioxins, PCBs, Other POPs

Pesticides in Food, Drugs in body fluid


EPA 1613, 8270A, 8290

# Magnetic Sector HRMS: Thermo Scientific DFS


- **Benchtop-like** operation (autotune)
- **Most compact** instrument in its class
- **Toroidal ESA** for highest **sensitivity**
- **Low power** consumption
- **Unique Dual GC** configuration
- **Fit-for-purpose** data evaluation software




# Magnetic Sector HRMS: Ion Optics of the DFS



# Magnetic Sector HRMS: Resolution Setting



**R 1000**




**R 10,000**

DFS: Tune peak of FC43 calibrant, m/z 414, screenshots at R 1,000 and R 10,000



# High Selectivity with 'High Resolution'

C12H4O2Cl4\*1.00 + C12H3Cl5\*1.00: p(gss, sfp.40) Chrg 10...



## Mass Peaks at Different Resolution Settings:

**Red:** TCDD mass trace at R 1.000

**Green:** TCDD mass trace at R 10.000

**Blue:** another chlorinated compound mass trace

**Accurate masses are used for target compound analysis**

# Need for High Resolution GC-MS – Blood Sample


## DFS GC-HRMS

Quantitation Mass  
 $m/z$  319.8965

*For the native TCDD*

**Chemical noise  
from matrix  
interferences  
shows up  
at lower mass  
resolution**


**ONLY HRMS at  
R=10,000  
selects the  
dioxin peak.**



# Triple Quadrupole TSQ8000 Evo for POPs analysis





**“Monitor a Transition”**  
from *Precursor* in Q1 to *Product Ion* in Q3



# Buffalo Milk Sample, TCDD Mass Traces


GC-MS/MS (5 µl injection PTV solv. split ) → GC-HRMS 1/10 concentration



Please note: Peak areas are not comparable between GC-MSMS and GC-HRMS

# Buffalo Milk Sample, PCDF

GC-MS/MS (5 µl injection PTV solv. split ) → GC-HRMS 1/10 concentration



Please note: Peak areas are not comparable between GC-MSMS and GC-HRMS


# TSQ 8000 Evo- TCDD in Buffalo Milk Samples

Blank(GC-MS/MS)


Buffalo Milk (GC-MS/MS)

Buffalo Milk (GC-HRMS)


2,3,7,8-TCDD




Labeled ISTD



0,17 pg/g fat



# Limit of Quantitation (LOQc) for GC-MS/MS




**Practically determined LOQc for PCDD/Fs**

**Serially diluted CS3 standards from CS3/100-CS3/1000**

**Ion ratio stability within  $\pm 15\%$  tolerance**


# TSQ 8000 Evo GC-MS/MS vs. GC-HRMS Analysis of CRM





# Sample preparation proposals


Courtesy: G Brambilla, ISS Rome



# TargetQuan - Quantitation Software Suite

## Features for dioxin applications:

- Toxicity equivalents (TEQ) - According to NATO and WHO definition, including
  - **Lower bound calculation:** requires using zero for the contribution of each non-quantified congener
  - **Medium bound calculation:** requires using half of the limit of quantification
  - **Upper boundary calculation:** requires using the limit of quantification
- User definable summation
  - Calculated amounts or TEQs
  - Reporting the sum TEQ values
- Isotopic ratio confirmation
  - One quantitation mass and up to two masses based on abundance
- EPA 1613 Rev.B compliant
  - Allows quantification based on average response of selected compound
  - Retention time correction



Entry	Compound Name	Entry Name	Entry Time	QM Retention Time	QM Area	RMF Area	RMF Ratio
1	12376-PeCDF	12376-PeCDF	27.52	12376-PeCDF	27.52	12376-PeCDF	1.0000
2	12376-PeCDF	12376-PeCDF	27.52	12376-PeCDF	27.52	12376-PeCDF	1.0000
3	12376-PeCDF	12376-PeCDF	27.52	12376-PeCDF	27.52	12376-PeCDF	1.0000
4	12376-PeCDF	12376-PeCDF	27.52	12376-PeCDF	27.52	12376-PeCDF	1.0000
5	12376-PeCDF	12376-PeCDF	27.52	12376-PeCDF	27.52	12376-PeCDF	1.0000
6	12376-PeCDF	12376-PeCDF	27.52	12376-PeCDF	27.52	12376-PeCDF	1.0000
7	12376-PeCDF	12376-PeCDF	27.52	12376-PeCDF	27.52	12376-PeCDF	1.0000
8	12376-PeCDF	12376-PeCDF	27.52	12376-PeCDF	27.52	12376-PeCDF	1.0000
9	12376-PeCDF	12376-PeCDF	27.52	12376-PeCDF	27.52	12376-PeCDF	1.0000
10	12376-PeCDF	12376-PeCDF	27.52	12376-PeCDF	27.52	12376-PeCDF	1.0000
11	12376-PeCDF	12376-PeCDF	27.52	12376-PeCDF	27.52	12376-PeCDF	1.0000
12	12376-PeCDF	12376-PeCDF	27.52	12376-PeCDF	27.52	12376-PeCDF	1.0000
13	12376-PeCDF	12376-PeCDF	27.52	12376-PeCDF	27.52	12376-PeCDF	1.0000
14	12376-PeCDF	12376-PeCDF	27.52	12376-PeCDF	27.52	12376-PeCDF	1.0000
15	12376-PeCDF	12376-PeCDF	27.52	12376-PeCDF	27.52	12376-PeCDF	1.0000
16	12376-PeCDF	12376-PeCDF	27.52	12376-PeCDF	27.52	12376-PeCDF	1.0000
17	12376-PeCDF	12376-PeCDF	27.52	12376-PeCDF	27.52	12376-PeCDF	1.0000
18	12376-PeCDF	12376-PeCDF	27.52	12376-PeCDF	27.52	12376-PeCDF	1.0000
19	12376-PeCDF	12376-PeCDF	27.52	12376-PeCDF	27.52	12376-PeCDF	1.0000
20	12376-PeCDF	12376-PeCDF	27.52	12376-PeCDF	27.52	12376-PeCDF	1.0000

# Applications and Methods

Application		Norm	Thermo Scientific TSQ 8000 Evo GC-MS/MS	Thermo Scientific DFS GC-HRMS
Food	EU Regulation	EU Regulatory Feed Control (at ML)	Approved	Approved
Food	EU Methods	EN 16215	Not in method	Approved
Food	EURL Recommendation	Background food studies (<1/5th EU ML)	Not recommended	Recommended
Clinical	EURL Recommendation	Human studies at trace levels	Not recommended	Recommended
Environmental	US Method	US EPA 8280A	Approved	Not in method
Environmental	EU Method	EN 1948	Not in method	Approved
Environmental	US Method	US EPA 1613 B for strict EPA compliance	Not in method	Approved
Environmental	US Method	US EPA Method 23	Not in method	Approved
Environmental	US Method	US EPA Method 8290	Not in method	Approved
Environmental	Japan Method	JIS K0311	Not in method	Approved
Environmental	Japan Method	JIS K0312	Not in method	Approved

- Any laboratory following EU approach will now be able to use GC-MS/MS for ML compliance control
- Laboratories strictly following US/Japan methods and regulation will continue to use GC-HRMS

## Conclusions - Why MS/MS and HRMS for POPs?

- **GC-MS/MS** – Versatile and Easy to use
  - Versatile instrumentation – no need to be dedicated
  - High sensitivity and selectivity for matrix samples
  - Provides TEQ results for PCDD/Fs and DL-PCBs
  - Fast - with high throughput capability, automated runs
- **GC-HRMS** – Unequivocal confirmation
  - Compliant with international regulations
  - Highest sensitivity and precision for ultra-low level analyses
  - Highly productive routine method – high throughput operation
- **GC-MS/MS + GC-HRMS** – Work seamlessly together
  - Use integrated sample prep workflow for screening and confirmation
  - High precision data using labelled internal standards
  - Common TargetQuan software platform



# Supplemental

**ThermoFisher SCIENTIFIC** Search All [dropdown] [search icon] [Contact Us] | Sign In [dropdown] Quick Order [shopping cart icon]

Recently Viewed Items [dropdown]

Applications & Techniques Shop All Products Services & Support About Us [cloud icon] Cloud

Home > Industrial & Applied Science > Environmental > Water Analysis > Wastewater Analysis > Instruments for Persistent Organic Pollutant (POPs) Analysis

## Instruments for Persistent Organic Pollutant (POPs) Analysis

← Wastewater Analysis

### Instruments for Persistent Organic Pollutant (POPs) Analysis

Automated Discrete Analyzers

Colorimeters



Thermo Fisher Scientific has a strong commitment to supporting laboratories charged with the task of researching and monitoring POPs. Our commitment ensures high productivity, added value solutions for some of the most challenging POPs determinations.

### Resources

- Power Plant Water Analysis Products
- Drinking Water Quality Analysis Process Products
- ROSS pH Electrodes Brochure
- SampleManager LIMS™ product brochure
- POPs Center of Excellence

### Support

Unity Lab Services—Access to our expert service team, information on support plans, parts and consumables resources, and other instrument services. [external link icon]

thermo scientific

applied biosystems

invitrogen

fisher scientific

unity lab services

# Sequence/Samples

- Two sequences were analysed on the system
  1. PCDD/Fs consisting of:-
    - CSL, CS1→CS6 EN:1948 standards in duplicate
    - Samples (fish and egg) bracketed by blank and CSL injections
    - Total sequence length 103 injection over ~4 days
  2. dl-PCBs consisting of: -
    - CS1→CS6 WP-CVS standards in duplicate
    - Samples (fish and egg) bracketed by CS4/100 diluted standards (up to 42 injections)
    - Serial dilutions of WP-CS4 standard (/100, /200, /400) to establish LOQ
    - Total sequence length 72 injections

# GC and MS Conditions

## TRACE 1310 GC Parameters

Injection Volume ( $\mu\text{L}$ ):	3
Liner:	LinerGOLD™ double taper (P/N: 453A1345-UI)
Inlet ( $^{\circ}\text{C}$ ):	280
Carrier Gas, (mL/min):	He, 1.2
Inlet Mode,	Splitless (split flow 120mL/min after 2 min)

## Oven Temperature Program:

Temperature 1 ( $^{\circ}\text{C}$ ):	170
Hold Time (min):	0
Temperature 2 ( $^{\circ}\text{C}$ ):	250
Rate ( $^{\circ}\text{C}/\text{min}$ ):	15.4
Hold Time (min):	0
Temperature 3 ( $^{\circ}\text{C}$ ):	285
Rate ( $^{\circ}\text{C}/\text{min}$ ):	2.5
Hold Time (min):	0
Temperature 3 ( $^{\circ}\text{C}$ ):	320
Rate ( $^{\circ}\text{C}/\text{min}$ ):	10
Hold Time (min):	15
Total Run Time (min):	39.7

## TSQ 8000 Evo Mass Spectrometer Parameters

Transfer Line ( $^{\circ}\text{C}$ ):	300
Ionization Type:	EI
Ion Source( $^{\circ}\text{C}$ ):	350
Electron Energy (eV):	40
Acquisition Mode:	Timed SRM with Dwell Time Prioritization
Acquired Masses (Da):	Optimised
Collision Energy (V):	Optimised



# Calibration Data PCDD/Fs

## EN:1948 CSL, CS1-CS6 PCDD/Fs

Congener (native)	R <sup>2</sup>	Mean RF	Std.Dev	RSD(%)
2378-TCDF	0.9999	0.999	0.029	2.9%
2378-TCDD	0.9998	1.131	0.048	4.3%
12378-PeCDF	0.9996	1.046	0.020	1.9%
23478-PeCDF	0.9999	1.086	0.026	2.4%
12378-PeCDD	0.9996	1.121	0.043	3.9%
123478-HxCDF	0.9998	1.153	0.025	2.1%
123678-HxCDF	1.0000	1.152	0.034	3.0%
234678-HxCDF	0.9997	1.117	0.025	2.2%
123478-HxCDD	0.9996	1.175	0.057	4.9%
123678-HxCDD	0.9995	1.087	0.029	2.7%
123789-HxCDD	0.9986	1.062	0.039	3.7%
123789-HxCDF	0.9990	1.058	0.035	3.3%
1234678-HpCDF	0.9999	1.078	0.022	2.1%
1234678-HpCDD	0.9998	1.070	0.034	3.2%
1234789-HpCDF	0.9999	1.112	0.023	2.1%
OCDD	0.9999	1.082	0.035	3.2%
OCDF	0.9998	1.518	0.038	2.5%

	1948CSL (pg/μl)	1948CS1 (pg/μl)	1948CS2 (pg/μl)	1948CS3 (pg/μl)	1948CS4 (pg/μl)	1948CS5 (pg/μl)	1948CS6 (pg/μl)
<b>NATIVE PCDDs &amp; PCDFs</b>							
2,3,7,8-Tetrachlorodibenzo-p-dioxin	0.04	0.2	0.8	4	16	80	320
1,2,3,7,8-Pentachlorodibenzo-p-dioxin	0.08	0.4	1.6	8	32	160	640
1,2,3,4,7,8-Hexachlorodibenzo-p-dioxin	0.08	0.4	1.6	8	32	160	640
1,2,3,6,7,8-Hexachlorodibenzo-p-dioxin	0.08	0.4	1.6	8	32	160	640
1,2,3,7,8,9-Hexachlorodibenzo-p-dioxin	0.08	0.4	1.6	8	32	160	640
1,2,3,4,6,7,8-Heptachlorodibenzo-p-dioxin	0.16	0.8	3.2	16	64	320	1280
Octachlorodibenzo-p-dioxin	0.16	0.8	3.2	16	64	320	1280
<b>SAMPLING STANDARDS</b>							
1,2,3,7,8-Pentachloro[ <sup>14</sup> C] <sub>12</sub> dibenzofuran	16	16	16	16	16	16	16
1,2,3,7,8,9-Hexachloro[ <sup>14</sup> C] <sub>12</sub> dibenzofuran	16	16	16	16	16	16	16
1,2,3,4,7,8,9-Heptachloro[ <sup>14</sup> C] <sub>12</sub> dibenzofuran	32	32	32	32	32	32	32
<b>EXTRACTION STANDARDS</b>							
2,3,7,8-Tetrachloro[ <sup>14</sup> C] <sub>12</sub> dibenzo-p-dioxin	16	16	16	16	16	16	16
1,2,3,7,8-Pentachloro[ <sup>14</sup> C] <sub>12</sub> dibenzo-p-dioxin	16	16	16	16	16	16	16
1,2,3,4,7,8-Hexachloro[ <sup>14</sup> C] <sub>12</sub> dibenzo-p-dioxin	16	16	16	16	16	16	16
1,2,3,6,7,8-Hexachloro[ <sup>14</sup> C] <sub>12</sub> dibenzo-p-dioxin	16	16	16	16	16	16	16
1,2,3,4,6,7,8-Heptachloro[ <sup>14</sup> C] <sub>12</sub> dibenzo-p-dioxin	32	32	32	32	32	32	32
Octachloro[ <sup>14</sup> C] <sub>12</sub> dibenzo-p-dioxin	32	32	32	32	32	32	32
<b>SYRINGE STANDARDS</b>							
1,2,3,4-Tetrachloro[ <sup>14</sup> C] <sub>12</sub> dibenzo-p-dioxin	16	16	16	16	16	16	16
1,2,3,7,8,9-Hexachloro[ <sup>14</sup> C] <sub>12</sub> dibenzo-p-dioxin	16	16	16	16	16	16	16


# Calibration Data dl-PCBs

## WP-CVS (WP-CS1 to WP-CS6) dl-PCBs

Congener (native)	R <sup>2</sup>	Mean RF	Std.Dev	RSD(%)
PCB -81- tetrachlorobiphenyl	0.9998	0.998	0.025	2.5%
PCB -77-tetrachlorobiphenyl	0.9999	1.032	0.027	2.6%
PCB -123- pentachlorobiphenyl	1.0000	0.937	0.011	1.2%
PCB -118 pentachlorobiphenyl	1.0000	0.993	0.016	1.6%
PCB -114 pentachlorobiphenyl	1.0000	1.037	0.027	2.6%
PCB -105 pentachlorobiphenyl	1.0000	0.956	0.020	2.1%
PCB - 126-pentachlorobiphenyl	0.9997	0.999	0.049	4.9%
PCB -167- hexachlorobiphenyl	0.9999	1.052	0.020	1.9%
PCB -156- hexachlorobiphenyl	0.9999	1.064	0.027	2.6%
PCB -157- hexachlorobiphenyl	0.9998	1.020	0.022	2.2%
PCB -169- hexachlorobiphenyl	0.9998	1.014	0.025	2.4%
PCB -189- heptachlorobiphenyl	1.0000	1.088	0.019	1.7%

NATIVE PCBs	IUPAC	WP-CS1 (ng/ml)	WP-CS2 (ng/ml)	WP-CS3 (ng/ml)	WP-CS4 (ng/ml)	WP-CS5 (ng/ml)	WP-CS6 (ng/ml)	WP-CS7 (ng/ml)
3,3',4,4'-Tetrachlorobiphenyl	77	0.1	0.5	2.0	10	40	200	800
3,4,4',5-Tetrachlorobiphenyl	81	0.1	0.5	2.0	10	40	200	800
2,3,3',4,4'-Pentachlorobiphenyl	105	0.1	0.5	2.0	10	40	200	800
2,3,4,4',5-Pentachlorobiphenyl	114	0.1	0.5	2.0	10	40	200	800
2,3',4,4',5-Pentachlorobiphenyl	118	0.1	0.5	2.0	10	40	200	800
2',3,4,4',5-Pentachlorobiphenyl	123	0.1	0.5	2.0	10	40	200	800
3,3',4,4',5-Pentachlorobiphenyl	126	0.1	0.5	2.0	10	40	200	800
2,3,3',4,4',5-Hexachlorobiphenyl	156	0.1	0.5	2.0	10	40	200	800
2,3,3',4,4',5'-Hexachlorobiphenyl	157	0.1	0.5	2.0	10	40	200	800
2,3',4,4',5,5'-Hexachlorobiphenyl	167	0.1	0.5	2.0	10	40	200	800
3,3',4,4',5,5'-Hexachlorobiphenyl	169	0.1	0.5	2.0	10	40	200	800
2,3,3',4,4',5,5'-Heptachlorobiphenyl	189	0.1	0.5	2.0	10	40	200	800
<b>MASS-LABELLED PCBs</b>								
3,3',4,4'-Tetrachloro[ <sup>13</sup> C <sub>12</sub> ]biphenyl	77L	50	50	50	50	50	50	50
3,4,4',5-Tetrachloro[ <sup>13</sup> C <sub>12</sub> ]biphenyl	81L	50	50	50	50	50	50	50
2,3,3',4,4'-Pentachloro[ <sup>13</sup> C <sub>12</sub> ]biphenyl	105L	50	50	50	50	50	50	50
2,3,4,4',5-Pentachloro[ <sup>13</sup> C <sub>12</sub> ]biphenyl	114L	50	50	50	50	50	50	50
2,3',4,4',5-Pentachloro[ <sup>13</sup> C <sub>12</sub> ]biphenyl	118L	50	50	50	50	50	50	50
2',3,4,4',5-Pentachloro[ <sup>13</sup> C <sub>12</sub> ]biphenyl	123L	50	50	50	50	50	50	50
3,3',4,4',5-Pentachloro[ <sup>13</sup> C <sub>12</sub> ]biphenyl	126L	50	50	50	50	50	50	50
2,3,3',4,4',5-Hexachloro[ <sup>13</sup> C <sub>12</sub> ]biphenyl	156L	50	50	50	50	50	50	50
2,3,3',4,4',5'-Hexachloro[ <sup>13</sup> C <sub>12</sub> ]biphenyl	157L	50	50	50	50	50	50	50
2,3',4,4',5,5'-Hexachloro[ <sup>13</sup> C <sub>12</sub> ]biphenyl	167L	50	50	50	50	50	50	50
3,3',4,4',5,5'-Hexachloro[ <sup>13</sup> C <sub>12</sub> ]biphenyl	169L	50	50	50	50	50	50	50
2,3,3',4,4',5,5'-Heptachloro[ <sup>13</sup> C <sub>12</sub> ]biphenyl	189L	50	50	50	50	50	50	50
<b>INTERNAL STANDARDS: MASS-LABELLED PCBs</b>								
2,3',4',5-Tetrachloro[ <sup>13</sup> C <sub>12</sub> ]biphenyl	70L	50	50	50	50	50	50	50
2,3,3',5,5'-Pentachloro[ <sup>13</sup> C <sub>12</sub> ]biphenyl	111L	50	50	50	50	50	50	50
2,2',3,4,4',5'-Hexachloro[ <sup>13</sup> C <sub>12</sub> ]biphenyl	138L	50	50	50	50	50	50	50
2,2',3,3',4,4',5-Heptachloro[ <sup>13</sup> C <sub>12</sub> ]biphenyl	170L	50	50	50	50	50	50	50

# Ion Ratio Stability



**Ion ratio stability for all 17 PCDD/F congeners over 103 consecutive sample/standard injections [n=16 EN:1948 CSL (40fg/ $\mu$ L 2,3,7,8-TCDD)]**

# Ion Ratio Stability

Congener	Ion ratio		Deviation	StDev	RSD(%)
	Theoretical	Calculated Mean			
2378-TCDF	0.96	0.93	-3%	0.063	6.7%
2378-TCDD	0.96	0.95	-2%	0.081	8.5%
12378-PeCDF	0.78	0.82	5%	0.047	5.8%
23478-PeCDF	0.78	0.81	3%	0.054	6.6%
12378-PeCDD	0.78	0.80	2%	0.067	8.4%
123478-HxCDF	0.63	0.62	-2%	0.039	6.3%
123678-HxCDF	0.63	0.63	0%	0.052	8.2%
234678-HxCDF	0.63	0.65	4%	0.039	6.0%
123478-HxCDD	0.64	0.63	-1%	0.054	8.6%
123678-HxCDD	0.63	0.64	1%	0.045	7.0%
123789-HxCDD	0.63	0.63	0%	0.047	7.5%
123789-HxCDF	0.63	0.65	3%	0.041	6.3%
1234678-HpCDF	0.80	0.78	-3%	0.056	7.3%
1234678-HpCDD	0.80	0.79	-1%	0.033	4.2%
1234789-HpCDF	0.80	0.77	-3%	0.034	4.3%
OCDD	0.96	0.94	-2%	0.079	8.4%
OCDF	0.96	0.97	1%	0.068	7.0%


**Ion ratio stability over 103 consecutive sample/standard injections (n=16), CSL ran after every 4 matrix sample injections.**

# Calculated Amount

Congener	Concentration (pg/ $\mu$ L)		Deviation	StDev	RSD(%)
	Nominal	Calculated Mean			
2378-TCDF	0.040	0.043	9%	0.003	7.1%
2378-TCDD	0.040	0.039	-2%	0.003	7.2%
12378-PeCDF	0.080	0.080	0%	0.003	3.8%
23478-PeCDF	0.080	0.079	-1%	0.005	5.9%
12378-PeCDD	0.080	0.077	-4%	0.006	7.8%
123478-HxCDF	0.080	0.084	5%	0.004	5.1%
123678-HxCDF	0.080	0.083	4%	0.004	4.7%
234678-HxCDF	0.080	0.079	-1%	0.005	5.8%
123478-HxCDD	0.080	0.078	-2%	0.007	8.4%
123678-HxCDD	0.080	0.083	4%	0.004	5.0%
123789-HxCDD	0.080	0.082	2%	0.007	8.0%
123789-HxCDF	0.080	0.081	2%	0.008	9.2%
1234678-HpCDF	0.160	0.171	7%	0.009	5.1%
1234678-HpCDD	0.160	0.178	11%	0.007	4.1%
1234789-HpCDF	0.160	0.164	3%	0.010	6.4%
OCDD	0.160	0.176	10%	0.014	7.8%
OCDF	0.160	0.176	10%	0.010	5.5%

Calculated concentration over 103 consecutive sample/standard injections (n=16), CSL ran after every 4 matrix sample injections.

# Limit of Confirmation (LOQc)




1948CS3 (pg/μl)	
<b>NATIVE PCDDs &amp; PCDFs</b>	
2,3,7,8-Tetrachlorodibenzo-p-dioxin	4
1,2,3,7,8-Pentachlorodibenzo-p-dioxin	8
1,2,3,4,7,8-Hexachlorodibenzo-p-dioxin	8
1,2,3,6,7,8-Hexachlorodibenzo-p-dioxin	8
1,2,3,7,8,9-Hexachlorodibenzo-p-dioxin	8
1,2,3,4,6,7,8-Heptachlorodibenzo-p-dioxin	16
Octachlorodibenzo-p-dioxin	16
2,3,7,8-Tetrachlorodibenzofuran	4
1,2,3,7,8-Pentachlorodibenzofuran	8
2,3,4,7,8-Pentachlorodibenzofuran	8
1,2,3,4,7,8-Hexachlorodibenzofuran	8
1,2,3,6,7,8-Hexachlorodibenzofuran	8
1,2,3,7,8,9-Hexachlorodibenzofuran	8
2,3,4,6,7,8-Hexachlorodibenzofuran	8
1,2,3,4,6,7,8-Heptachlorodibenzofuran	16
1,2,3,4,7,8,9-Heptachlorodibenzofuran	16
Octachlorodibenzofuran	16
<b>SAMPLING STANDARDS</b>	
1,2,3,7,8-Pentachloro[ <sup>13</sup> C <sub>12</sub> ]dibenzofuran	16
1,2,3,7,8,9-Hexachloro[ <sup>13</sup> C <sub>12</sub> ]dibenzofuran	16
1,2,3,4,7,8,9-Heptachloro[ <sup>13</sup> C <sub>12</sub> ]dibenzofuran	32
<b>EXTRACTION STANDARDS</b>	
2,3,7,8-Tetrachloro[ <sup>13</sup> C <sub>12</sub> ]dibenzo-p-dioxin	16
1,2,3,7,8-Pentachloro[ <sup>13</sup> C <sub>12</sub> ]dibenzo-p-dioxin	16
1,2,3,4,7,8-Hexachloro[ <sup>13</sup> C <sub>12</sub> ]dibenzo-p-dioxin	16
1,2,3,6,7,8-Hexachloro[ <sup>13</sup> C <sub>12</sub> ]dibenzo-p-dioxin	16
1,2,3,4,6,7,8-Heptachloro[ <sup>13</sup> C <sub>12</sub> ]dibenzo-p-dioxin	32
Octachloro[ <sup>13</sup> C <sub>12</sub> ]dibenzo-p-dioxin	32
2,3,7,8-Tetrachloro[ <sup>13</sup> C <sub>12</sub> ]dibenzofuran	16
2,3,4,7,8-Pentachloro[ <sup>13</sup> C <sub>12</sub> ]dibenzofuran	16
1,2,3,4,7,8-Hexachloro[ <sup>13</sup> C <sub>12</sub> ]dibenzofuran	16
1,2,3,6,7,8-Hexachloro[ <sup>13</sup> C <sub>12</sub> ]dibenzofuran	16
2,3,4,6,7,8-Hexachloro[ <sup>13</sup> C <sub>12</sub> ]dibenzofuran	16
1,2,3,4,6,7,8-Heptachloro[ <sup>13</sup> C <sub>12</sub> ]dibenzofuran	32
Octachloro[ <sup>13</sup> C <sub>12</sub> ]dibenzofuran	32
<b>SYRINGE STANDARDS</b>	
1,2,3,4-Tetrachloro[ <sup>13</sup> C <sub>12</sub> ]dibenzo-p-dioxin	16
1,2,3,7,8,9-Hexachloro[ <sup>13</sup> C <sub>12</sub> ]dibenzo-p-dioxin	16

Figure 2. Practically determined LOQc for PCDD/Fs

Serially diluted CS3 standards from CS3/100-CS3/1000




Ion ratio stability within ±15% tolerance

# Chromatography – Lowest Calibration Standard (EN:1948 CSL)



Natives only

# Sample – Selected Chromatography And Calculations



123789-HxCDF – “below” LOQc

12378-PeCDD – 0.162 pg/g WHO-TEQ


2378-TCDD – 0.193 pg/g WHO-TEQ

Peak Name	Retention Time min	Amount pg/ $\mu$ L	LOQc pg/ $\mu$ L	WHO-TEF (2005)	WHO-TEQ(2005) pg/g	Peak Confirmation
2378-TCDF	18.642	12.8690	0.02	0.1	0.512	Valid
2378-TCDD	19.103	0.4864	0.02	1	0.193	Valid
12378-PeCDF	21.977	0.7828	0.048	0.03	0.009	Valid
23478-PeCDF	22.956	2.1316	0.04	0.3	0.254	Valid
12378-PeCDD	23.117	0.4070	0.072	1	0.162	Valid
123478-HxCDF	25.356	0.1276	0.04	0.1	0.005	Valid
123678-HxCDF	25.456	0.1242	0.048	0.1	0.005	Valid
234678-HxCDF	26.056	0.1160	0.048	0.1	0.005	Valid
123478-HxCDD	26.118	0.0938	0.064	0.1	0.004	Valid
123678-HxCDD	26.207	0.2346	0.064	0.1	0.009	Valid
123789-HxCDD	26.507	0.0825	0.072	0.1	0.003	Valid
123789-HxCDF	26.878	0.0301	0.064	0.1	0.001	nd<LOQc
1234678-HpCDF	28.396	0.1036	0.112	0.01	0.000	nd<LOQc
1234678-HpCDD	29.747	0.4449	0.112	0.01	0.001	Valid
1234789-HpCDF	30.475	0.0453	0.096	0.01	0.000	nd<LOQc
OCDD	34.197	0.4003	0.112	0.0003	0.000	Valid
OCDF	34.457	0.0417	0.096	0.0003	0.000	nd<LOQc
<b>Sum</b>					<b>1.163</b>	<b>pg/g</b>

Fish sample WHO-TEQ (2005) – 78.6g Fresh weight – HR-GCMS 1.29pg/g




# Sample Chromatography And Calibration – PCDD/Fs



Fish sample WHO-TEQ (2005) 1.16pg/g – 78.6g Fresh weight

# Standard Chromatography And Calibration – dl-PCBs (WP-CS1)



# Dioxin Analysis Legislation - Approaches to Dioxin Analysis

- Europe
  - Performance based approach.
  - Regulated technology and defined performance and quality criteria.
  - Methods are provided as a guide and can be modified.
  - Levels in food are regulated with strict maximum limits.
  - Emission levels are regulated.
- North America & Japan
  - Method based approach
  - Regulatory agencies (such as USEPA) provide complete methods to be followed.
  - Deviations to the method are possible. But must be documented closely.
    - Significant deviations from methods can result in the analysis being challenged or becoming unsellable in market.
  - Levels in food are not regulated.
  - Emission guidelines are given.
- ROW adopt variety of approaches,
  - dependent on regulatory and/or business requirement

## Background Info: Confirmation vs Screening

- **Confirmatory Analysis**

- Unequivocal identification and quantification
- Results can be used in prosecution



- **Screening Analysis**

- Indication of levels only
- Must be followed with confirmatory analysis



# What Has Been Happening In EU Dioxin Analysis?

## ANALYTICAL CRITERIA FOR USE OF MS/MS FOR DETERMINATION OF DIOXINS AND DIOXIN-LIKE PCBS IN FEED AND FOOD

Kotz A<sup>1</sup>, Malisch R<sup>1\*</sup>, Focant J<sup>2</sup>, Eppe G<sup>2</sup>, Cederberg TL<sup>3</sup>, Rantakokko P<sup>4</sup>, Fürst P<sup>5</sup>, Bernsmann T<sup>5</sup>, Leondiadis L<sup>6</sup>, Lovász C<sup>7</sup>, Scortichini G<sup>8</sup>, Diletti G<sup>8</sup>, di Domenico A<sup>9</sup>, Ingelido AM<sup>9</sup>, Traag W<sup>10</sup>, Smith F<sup>11</sup>, Fernandes A<sup>11</sup>

<sup>1</sup>European Union Reference Laboratory (EU-RL) for Dioxins and PCBs in Feed and Food, Freiburg, Germany; <sup>2</sup>CART, University of Liege, Belgium; <sup>3</sup>National Food Institute, Technical University of Denmark, Søborg, Denmark; <sup>4</sup>National Institute for Health and Welfare (THL), Kuopio, Finland; <sup>5</sup>Chemisches und Veterinäruntersuchungsamt Münsterland-Emscher-Lippe (CVUA-MEL), Münster, Germany; <sup>6</sup>NCSR Demokritos, Athens, Greece; <sup>7</sup>Central Agricultural Office, Budapest, Hungary; <sup>8</sup>Istituto Zooprofilattico Sperimentale dell'Abruzzo e del Molise "G. Caporale, Teramo, Italy; <sup>9</sup>Istituto Superiore di Sanità (ISS), Roma, Italy; <sup>10</sup>RIKILT, Wageningen, The Netherlands; <sup>11</sup>The Food and Environment Research Agency (FERA), York, United Kingdom

<http://www.dioxin20xx.org/pdfs/2012/1041.pdf>

- 2009–2012 GC-MS/MS systems were evaluated from several vendors
- Conclusion:
  - “GC-MS/MS systems with sufficient sensitivity” can be “used as confirmatory methods”
  - “However, for the determination of low background levels the application of GC-HRMS is required” e.g. <20% of Maximum Level (ML)

# Does MS/MS Agree with HRMS?




Figure 11. Comparative results (upperbound concentration values) for 50 food and feed samples analyzed by GC/HRMS and GC/MS/MS.

- At EU ML results are comparable.
- Below EU ML results are not comparable.


# Does MS/MS Agree with HRMS?

## Comparison of GC-MS/MS Results of laboratories (3)

PCDD/Fs

- Chemisches und Veterinäruntersuchungsamt (CVUA) Freiburg (EU-RL for Dioxins and PCBs), GC-MS/MS

Comparable results for MS/MS and HRMS for WHO-TEQ, deviation of MS/MS mostly below 20 %



14 / GC-MS/MS

**EURL**  
European Union Reference Laboratory  
Dioxins and PCBs

State Institute for Chemical and  
Veterinary Analysis of Food  
CVUA Freiburg



- At EU ML results are comparable.
- Below EU ML results are not comparable.

# What's New In EU Regulation For Dioxins?

## Past

[Commission Regulation \(EU\) No 252/2012 of 21 March 2012](#)

- Specifies use of **GC-HRMS** for confirmatory dioxin analysis
- **GC-MS/MS** was allowed a screening technique.



## Present

[Commission Regulation \(EU\) No 589/2014 of 2 June 2014](#)

- Specifies use of **GC-HRMS** or **GC-MS/MS** for confirmatory dioxin analysis
- **GC-MS/MS** is “an appropriate confirmatory method for checking compliance with the maximum”, only.
- **GC-HRMS** remains the recommended technique for “determination of low background levels in food monitoring, following of time trends, exposure assessment of the population”.



# Compliance With New EU Dioxin Regulation

European Commission GC-MS/MS Confirmatory Performance Criteria	TSQ 8000 Evo GC-MS/MS Capabilities	Compliance Confirmed
Two specific precursor ions with two specific production ions	All recommended methods developed <b>as defined in criteria</b> Fully optimized by AutoSRM	<b>Yes</b>
Tolerance of ion ratios within $\pm 15\%$	<b>&lt; 10%</b> measured at EPA 1613 CSL level (n=14)	<b>Yes</b>
Resolution of each quadrupole equal to or better than unit mass resolution	All recommended methods developed <b>Q1 and Q3 at 0.7 Da</b>	<b>Yes</b>
The % RSD of the five (or more) Relative Response Factors (RRFs) for each unlabeled PCDD/PCDF and labelled internal standards must not exceed 20%	6 point curve EPA 1613 CSL-CS4 <b>2%</b>	<b>Yes</b>

# Applications and Methods

Application		Norm	Thermo Scientific TSQ 8000 Evo GC-MS/MS	Thermo Scientific DFS GC-HRMS
Food	EU Regulation	EU Regulatory Feed Control (at ML)	Approved	Approved
Food	EU Methods	EN 16215	Not in method	Approved
Food	EURL Recommendation	Background food studies (<1/5th EU ML)	Not recommended	Recommended
Clinical	EURL Recommendation	Human studies at trace levels	Not recommended	Recommended
Environmental	US Method	US EPA 8280A	Approved	Not in method
Environmental	EU Method	EN 1948	Not in method	Approved
Environmental	US Method	US EPA 1613 B for strict EPA compliance	Not in method	Approved
Environmental	US Method	US EPA Method 23	Not in method	Approved
Environmental	US Method	US EPA Method 8290	Not in method	Approved
Environmental	Japan Method	JIS K0311	Not in method	Approved
Environmental	Japan Method	JIS K0312	Not in method	Approved

- Any laboratory following EU approach will now be able to use GC-MS/MS for ML compliance control
- Laboratories strictly following US/Japan methods and regulation will continue to use GC-HRMS