


**Thermo Fisher Scientific
OES and XRF/XRD product range**
Sofia 2017
ACM2 Anniversary
Peter Bellant

Elemental Spectroscopy - Detection Limits and Dynamic Range




Legend

- Solids
- Liquids, gases

Arc/Spark Optical Emission Spectrometry

- In Optical Emission Spectrometry
 - We observe low energy electronic transitions
 - Typically several eV ($1 \text{ eV} \cong 23 \text{ Kcal mol}^{-1}$)
 - From VUV to NIR (100 - 900nm)



Why using Arc/Spark OES?

- Extremely fast and simple analysis of conducting solids
- Can analyze simultaneously > 40 elements from trace (ppm) to major (%) levels in <1 minute
 - ppm or sub-ppm limits of detections
 - Can analyze P, S, C, N, O at low levels
- Sample preparation is fast and simple
- Accurate, precise, stable, reliable
- Cheap analysis: low costs of ownership and maintenance
- High instrument availability
- Long life and robustness
- Well established (exists since 1934)

Use of Spark-OES

- Users of OES

- Primary metals producers
- Foundries
- Metals processors
- Metals recyclers
- Service laboratories
- Universities, research centers



- OES instruments main uses

- Process control
- Final product quality control
- Research
- Certificate of analysis
- Incoming materials control
- Metal sorting

- Prerequisites for Spark-OES analysis


- Solid conductive sample
- Representative sample (sample taking)
- Clean and flat sample surface (sample preparation)
- Reference materials for calibration (CRM and RM)
- Setting-up samples (SUS) for drift correction (to maintain accuracy)

Applications of Optical Emission Spectrometry

Elements and matrices analyzable by OES

ARC/SPARK OPTICAL EMISSION SPECTROMETRY (OES) PERIODIC CHART

ELEMENTS AND MATRICES ANALYZABLE BY OES




www.thermo.com/oes

14 routinely analyzed matrices

7 special matrices

5 potential matrices

Spark-OES instrument operational diagram



ARL 3460

Mid-range to high-end PMT-based instrument

- The market reference
 - > 6'000 units sold since launch
- Mono- and multi-base
- Up to 60 channels in standard
- Standard Hirep source
- Conventional acquisition electronics
- Key differentiators
 - Reputation, large installed base
 - Overall performance, including stability
 - Reliability, robustness, lifetime
 - Quality of the calibration
 - OXSAS analytical software with SPC, Total Materia (Key to Metals) included, measurement uncertainty, Total Materia (Key to Metals) option, etc.



ARL 4460

Ultimate performance PMT instrument

- For many OES users, the best instrument on the market
- >2'300 units sold – a real success!
- Best sensitivity, speed of analysis, precision, accuracy on the market
- Digital technologies
 - Excitation with Current Controlled Source (CCS)
 - Acquisition with Time Resolved Spectroscopy (TRS)
- Inclusions analysis with the Spark-DAT option
 - Probably the fastest method for inclusion analysis
 - Can be used to control inclusions during production



ARL iSpark 8860 – Experience and performance

ARL
iSPARK 8860
OES Metals Analyzer



- Typical high-performance single matrix PMT instrument
 - Customized up to 80 PMTs
 - Typically for demanding steel plants, pure metals producers, companies that need best possible analysis of all the elements at trace level
- Some benefits
 - High level ultra-pure metal and trace element analysis
 - High analysis speed
 - Advanced signal acquisition and processing methods for all the channels
 - Highly reliable, high throughput automated analysis
 - Advanced ultra-fast inclusion analysis options available
 - Easy and fast stand maintenance for increased instrument availability
 - Synoptic tool for fast visual trouble-shooting
 - Significant argon savings

ARL easySpark – Technical Highlights

- Innovative high safety open stand
- Multi grating / CCD spectrograph allowing continuous and wide elemental coverage including UV elements like Nitrogen
- Compact “intelliSource” CCS Source
- Robust, stable and low-frequency, easy maintenance stand
- Optimized argon stream design for improved analysis performances on difficult samples
- EasyOXAS software with user-friendly interface matching all metal producers needs





ThermoFisher
SCIENTIFIC

X-Ray Fluorescence X-ray Diffraction

Applications and products

Thermo Fisher Scientific
1024 Ecublens, Switzerland



The world leader in serving science

What is so great about WD X-ray Fluorescence ?

- Precise, non-destructive multi-element analysis of conducting and non-conducting samples:
 - Solids
 - Liquids
 - Loose powders
- Multi-element analysis: Be - F to U
- Inorganic and Organic materials
- High precision and highly reliable
- Wide dynamic range: sub-ppm to 100%
- Ten times better resolution compared to EDXRF (15eV at Mn $K\alpha$)
- Better performance on light elements compared to EDXRF (better stability and limits of detection)
- Powerful “**standard-less**” elemental analysis for totally unknown samples

WDXRF, EDXRF and XRD: Elemental and Phase analysis of a variety of materials

- *Cement and building materials*
- *Metals, Slags, ferro-alloys*
- *Petroleum, Polymers, Oils*
- *Ores and raw materials*
- *Chemicals/Pharmaceuticals*
- *Geology, geochemistry*
- *Environmental, air filters*
- *Food products*
- *Mining extraction*
- *Universities, central labs*
- *Thin films, magnetic media, paints*
- *Etc.*



Thermo Scientific XRF and XRD Product Portfolio: Strong and complementary technologies

EDXRF



ARL QUANT'X
Advanced EDXRF

WDXRF



ARL PERFORM'X
High Performance
sequential XRF

Integrated XRF and XRD



ARL 9900 Series
Integrated XRF-XRD

Powder XRD



Equinox 6000: High
performance Powder XRD



Portable Niton



ARL OPTIM'X: Surprising
performance in WDXRF



Equinox 100 & 1000
Benchtop XRD

XRF: Elemental analysis

**XRD: Analysis of
Structure-crystallography
Phase or compound**

Thermo Scientific EDXRF Spectrometer

ARL QUANT'X



- Most suitable for
 - Gemology
 - RoHS and WEEE compliance analysis
 - Forensics and Trace Analysis
 - Environmental analysis
 - Oils and gasoline
 - Polymers and plastics
 - Industrial minerals
 - Nutritional Supplements in Food
 - Aerosol Particulate Filters
 - Alternative fuels and backup in cement
 - Slags in steel industry

Thermo Scientific Wavelength Dispersive XRF systems

ARL OPTIM'X XRF

- Compact WDXRF instrument with a unique combination of sequential and simultaneous devices
- Amazing performance despite low power thanks to UCCO technology
- Choice of 50W or 200W power
- Excellent spectral resolution
 - Better accuracy and precision of analysis
- Markets
 - Raw materials
 - Slags, limestone, sand, feldspar, etc.
 - Ferro-alloys, aluminium or zinc ingots, etc.
 - Hot metal, irons
 - Food and Chemicals
 - Back-up in big cement and steel plants
 - Refineries (Sulfur Analyzer)



ARL Optim'X
with ARL SMS-Omega

ARL PERFORM'X series

- Highest performance sequential XRF spectrometers for qualitative, quantitative and standard-less analysis of almost any material
- Applications:
 - Central Laboratories, Research Labs, Materials Science
 - Geology centers, Petrochemicals, Polymers,
 - Environmental and process control



4200W

2500W

1500W



- **New X-ray tube**
 - 50 μ Be window
 - Extreme stability



**Versatile
Goniometer**



**Integrated Large XY
Sample changer**

- + Mapping/Spots down to 0.5mm
- + Small samples down to 5mm
- + Helium for liquid/loose powders

Thermo WDXRF systems

ARL 9900 Series X-ray Spectrometer



- Compact, modular XRF instrument for rapid elemental sim/seq analysis
 - Up to 32 fixed channels
 - Up to 24 fixed channels with one universal goniometer
 - Dual goniometer + 8 fixed channels
- Applications:
 - Iron and Steel
 - Aluminum
 - Copper
 - Cement
 - Mining and metal extraction

Comparison of XRF and OES techniques: Typical precision on various elements in Steel

		ARL 9900 XRF	ARL 4460 OES
Element	Conc.	in 30s ppm	ppm
Si	0.002	2.6	0.6
Si	0.5	8.2	22
Si	1	11.4	40
Si	5	25.0	150
Mn	0.001	2	0.5
Mn	0.1	3.2	5
Mn	1	8.4	30
Mn	5	18.4	130
Mn	10	26.6	250
P	0.01	1.4	1
P	0.3	5.5	20
S*	0.01	1.3	2.5
S*	0.05	2.6	10

*Low Mo samples

		ARL 9900 XRF	ARL 4460 OES
Element	Conc.	in 30s ppm	ppm
Cr	0.01	1.3	1
Cr	0.1	2.5	4
Cr	1	7.2	25
Cr	5	16.0	80
Cr	10	23	150
Al	0.001	1.7	0.3
Al	0.01	2.4	1
Al	0.1	5.5	9
Ni	0.001	2	0.8
Ni	0.1	4	3
Ni	1	10	25
Ni	10	32	250
Ni	30	110	750

ARL 9900 XRF – 3.6kW compared to ARL 4460 OES:
XRF is better at high levels – OES is better at low levels

Thermo Scientific Combined XRF-XRD systems

ARL 9900 Series X-ray Spectrometer




- Compact, modular XRF instrument for rapid elemental sim/seq analysis
- Specific phase analysis using **unique** patented integrated XRD system (e.g. Total Cement Analyzer - TCA)
- Applications for XRD system:
 - Cement
 - Free lime and clinker phases
 - Hot meal decarbonation
 - Limestone additions to cement
 - Clinker content in cement
 - Iron and Steel
 - Fe 2+ in sinters
 - Free lime in slags
 - DRI
 - Aluminum
 - Bath ratio of electrolytic bath
 - Copper mining and metal extraction
 - Various phases (compounds)

NEW offer on XRD: EQUINOX 100 & 1000

EQUINOX 100 & 1000

Bench-top X-ray Diffractometers

- Ultra-fast unique Equinox curved detector
 - measures all diffraction peaks simultaneously
 - also called PSD for position sensitive detector



- Unique simultaneous XRD detection



Environment



Geology, Cement, Mining



Research



Metallurgy



Environment, Geology



Pharmaceuticals



Education



Polymers, Biomaterials



Electronics



Photovoltaic