




Recent developments in ICP-MS - Introducing the Thermo Scientific iCAP TQ ICP-MS

Sofia 18.10.2017


Burkhard Stehl, Sales Manager Trace Elemental Analyzers ACO

Introducing our ICP-MS portfolio


Technology for all challenges




Thermo Scientific™ NEPTUNE Plus™ MC-ICP-MS
Multicollector ICP-MS



Thermo Scientific™ ELEMENT2/XR™ HR-ICP-MS
High Resolution ICP-MS



Thermo Scientific™ iCAP™ TQ ICP-MS
Triple quadrupole ICP-MS




Thermo Scientific™ iCAP™ RQ ICP-MS
Single quadrupole ICP-MS




Introducing our ICP-MS portfolio

Technology for all challenges




Thermo Scientific™ NEPTUNE Plus™ MC-ICP-MS
Multicollector ICP-MS



Thermo Scientific™ ELEMENT2/XR™ HR-ICP-MS
High Resolution ICP-MS



Thermo Scientific™ iCAP™ TQ ICP-MS
Triple quadrupole ICP-MS



Thermo Scientific™ iCAP™ RQ ICP-MS
Single quadrupole ICP-MS




All the Power, None of the Complexity

- ✓ **Advanced interference removal**
- ✓ **Robust design for routine analysis**
- ✓ **Integrated automation options**
- ✓ **Flexible for advanced applications**
- ✓ **Unique ease of use – Reaction Finder**


**Triple quadrupole accuracy with
single quadrupole ease of use**



iCAP TQ ICP-MS: How it works - on mass reaction mode




iCAP TQ ICP-MS: How it works - product ion reaction mode



All the Power, None of the Complexity

- Problem: when faced with measurement of a sample where interferences are expected, which is the best measurement mode?
- Solution: method development assistant – **intelligent Reaction Finder, iRF**
 - Software concept for intelligent selection of all 3 parameters
 - Just select the element for analysis and the software does the rest



Reaction Finder method development assistant

Without Reaction Finder

Select

- Select the Analytes to be measured

Select

- For each analyte, select the isotopes to be measured

Select

- Select the internal standard element

Select

- Select the Q1 Analyte

Select

- Select the CRC gas (None, He, H₂, O₂, NH₃)

Select

- Select the mode (KED, Single Quad Mode, Triple Quad Mode)

Select

- Select the Q3 Mass (On-mass/mass shift product ion)

Decide

- Are the suggested settings ok? If not, update them

Analyze

- Enter sample names and positions or import from LIMS and start the LabBook

With Reaction Finder

Select

- Select the Analytes to be measured

Select

- Select the internal standard element

Decide

- Are the suggested settings ok? If not, update them

Analyze

- Enter sample names and positions or import from LIMS and start the LabBook



Reaction Finder in Thermo Scientific™ Qtegra™ ISDS Software

Reaction Finder is a supplied applet that preselects optimised conditions for each target isotope in each available mode


For example for ^{31}P , the Reaction Finder database defines the following method parameters:

Analyte type	Analyte	Is default isotope	Reaction gas	Q1 mass (u)	Q3 analyte	Is default Q3 Analyte	Is default reaction
Isotope	31P	<input checked="" type="checkbox"/>	O ₂ (Oxygen)	30.9737634	31P	<input type="checkbox"/>	<input type="checkbox"/>
Isotope	31P	<input checked="" type="checkbox"/>	O ₂ (Oxygen)	30.9737634	31P.16O	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Isotope	31P	<input checked="" type="checkbox"/>	O ₂ (Oxygen)	30.9737634	31P.17O	<input type="checkbox"/>	<input type="checkbox"/>
Isotope	31P	<input checked="" type="checkbox"/>	O ₂ (Oxygen)	30.9737634	31P.18O	<input type="checkbox"/>	<input type="checkbox"/>
Isotope	31P	<input checked="" type="checkbox"/>	O ₂ (Oxygen)	30.9737634	31P.16O2	<input type="checkbox"/>	<input type="checkbox"/>
Isotope	31P	<input checked="" type="checkbox"/>	O ₂ (Oxygen)	30.9737634	31P.17O.16O	<input type="checkbox"/>	<input type="checkbox"/>
Isotope	31P	<input checked="" type="checkbox"/>	O ₂ (Oxygen)	30.9737634	31P.18O.16O	<input type="checkbox"/>	<input type="checkbox"/>
Isotope	31P	<input checked="" type="checkbox"/>	O ₂ (Oxygen)	30.9737634	31P.17O2	<input type="checkbox"/>	<input type="checkbox"/>
Isotope	31P	<input checked="" type="checkbox"/>	O ₂ (Oxygen)	30.9737634	31P.18O.17O	<input type="checkbox"/>	<input type="checkbox"/>
Isotope	31P	<input checked="" type="checkbox"/>	O ₂ (Oxygen)	30.9737634	31P.18O2	<input type="checkbox"/>	<input type="checkbox"/>
Isotope	31P	<input checked="" type="checkbox"/>	H ₂ (Hydrogen)	30.9737634	31P	<input type="checkbox"/>	<input type="checkbox"/>
Isotope	31P	<input checked="" type="checkbox"/>	H ₂ (Hydrogen)	30.9737634	31P.1H4	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Isotope	31P	<input checked="" type="checkbox"/>	None (No reaction gas)	30.9737634	31P	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Isotope	31P	<input checked="" type="checkbox"/>	He (Helium)	30.9737634	31P	<input checked="" type="checkbox"/>	<input type="checkbox"/>


None of the complexity, all of the flexibility:

- **Default reactions for all modes of iCAP TQ ICP-MS operation including collision/ reaction gases such as O₂, H₂, NH₃ and He**
- **Dedicated mass flow controller for each cell gas**

iCAP TQ ICP-MS – Feature summary



Intuitive quick-connect sample introduction components



Interface design

Bench level pop-out interface for easy ambidextrous access to the cones

and

the extraction lens for simplest possible routine maintenance

...without needing to break the vacuum



Ion focusing: the *RAPID* lens

Right

Angle

Positive

Ion

Deflection


90° ion focusing with
total ion deflection in 3
dimensions

and

Elimination of neutral
species

for

Highest signal to noise
ratio of any ICP-MS



The background of the slide features a large, semi-transparent globe of the Earth in the center, surrounded by several smaller, textured blue spheres of varying sizes. The overall color palette is dominated by shades of blue and teal, with a white horizontal band across the middle containing the logo and title, and a solid red horizontal band at the bottom.

ThermoFisher
SCIENTIFIC

Applications - iCAP TQ ICP-MS

The world leader in serving science

Main application areas for triple quadrupole ICP-MS

Meeting human health
and environmental challenges

Advancing development in
metals, materials and chemicals



- **Clinical Research and Toxicology**
- **Metallopharmaceuticals**
- **Environmental Analysis/Monitoring**
- **Food Safety**

- **Material Analysis**
- **Nanoparticle Characterization**
- **Metallurgy**
- **Energy Production**

iCAP TQ measurement modes

- SQ – mode – $\text{H}_2/\text{He}/\text{KED}/\text{O}_2/\text{NH}_3$
- TQ – mode – $\text{He}/\text{H}_2/\text{O}_2/\text{NH}_3$
- Product ion measurement (analyte ion is reactive and moved to a new product ion mass).
 - $^{32}\text{S} + ^{16}\text{O} \rightarrow ^{32}\text{S}^{16}\text{O}$
- On mass measurement (interfering ions are reactive and moved away from analyte ion).
 - ^{172}Yb Use NH_3 to remove $^{156}\text{Gd}^{16}\text{O}$



Spectral overlaps

- Avoidance of spectral overlaps on reaction product ions with O₂ cell gas
- SQ – ICP-MS – can use same reactive chemistry
- TQ ICP-MS – filtration of ions prior to QCell
- Spectral overlaps on product ions – can cause errors in results.


Q1	Q3	Overlaps		
		Ni	Cu	Zn
Ti	TiO			
⁴⁶ Ti	62 - ⁴⁶ Ti ¹⁶ O	⁶² Ni		
⁴⁷ Ti	63 - ⁴⁷ Ti ¹⁶ O		⁶³ Cu	
⁴⁸ Ti	64 - ⁴⁸ Ti ¹⁶ O			⁶⁴ Zn
⁴⁹ Ti	65 - ⁴⁹ Ti ¹⁶ O		⁶⁵ Cu	
⁵⁰ Ti	66 - ⁵⁰ Ti ¹⁶ O			⁶⁶ Zn

Example - Titanium


- Prepare 1ppb Ti solution spiked with 10ppb Ni, Cu and Zn
- Measure using SQ and TQ modes with oxygen in QCell
- Set up scan to view mass region 60-70
- Overlay spectra to compare data from SQ and TQ modes.



1ppb Ti – TQ and SQ modes (Oxygen)



1ppb Ti – TQ – mass shift oxygen



Arsenic and selenium in environmental samples



Low levels in samples



Multiple interferences:
 ArCl^+ , Ar_2^+ and REE^{++}


33
As
Arsenic
74.921595

34
Se
Selenium
78.971




As and Se analysis in the presence of REE's – the problem

Usual interferences on As and Se - Ar₂, ArCl - easy to remove using He KED, but if REE are present...



As and Se analysis in the presence of REE's: the iCAP TQ solution


- Control ions entering the collision cell using **Q1**
- Use O_2 to efficiently convert As and Se to AsO^+ and SeO^+ in **Q2** (i.e. the collision cell)
- REE⁺⁺ species don't react
- Selectively detect AsO^+ (at mass 91) and SeO^+ (at mass 94) free from REE⁺⁺ interference, using **Q3**



Type	^{75}As	Method to remove	^{78}Se	Method to remove
Polyatomic	$^{40}Ar^{35}Cl$	KED	$^{40}Ar^{38}Ar$	KED, H_2
	$^{40}Ca^{35}Cl$			
Isobaric	$^{150}Nd^{2+}$	O_2	$^{156}Gd^{2+}$	O_2
	$^{150}Sm^{2+}$		$^{156}Dy^{2+}$	

As and Se with REE present - results in different modes

Interference removal capability in each mode



- 1ppm Dy, Gd, Nd, Sm and Tb added
- Increased BECs observed for all SQ-modes due to unresolved doubly charged REE interferences
- Hydrogen is suitable for removing Ar based polyatomics, but is not capable of fully removing REE²⁺ interferences
- TQ-O₂ mode shows dramatically lower BEC values for both As and Se
- Accuracy assessed by analysis of AGV andesite reference material and a deep sea sediment
- Spike recovery tests also performed

How do we know iMS is effective?

- Let's look at an example using Se in the presence of high As

Sample	Added Signal for Hypothetical $^{75}\text{As}^{18}\text{O}$ (cps)	Calculated $^{77/80}\text{Se}$ Isotope Ratio	Measured $^{77/80}\text{Se}$ Isotope Ratio
0.5 ppm Se	+0	0.1398	0.1398
0.5 ppm Se, 50 ppb As	+1,260	0.1402	0.1394
0.5 ppm Se, 100 ppb As	+2,493	0.1410	0.1392
0.5 ppm Se, 250 ppb As	+7,298	0.1443	0.1393
0.5 ppm Se, 500 ppb As	+14,791	0.1497	0.1391

Measured isotope ratios are not corrected for mass bias

**Addition of As to the sample solution has no effect on the $^{77/80}\text{Se}$ ratio
→ iMS is working perfectly well!**

Determination of Ti in biological samples using ICP-MS

22

Ti

Titanium


47.867

Titanium based components used for orthopedic and dental implants.

Degradation of these implants releases Ti (and Co, Ni and Cr too) into the body


$^{48}\text{Ca}^+$, PO^+ , SO^+ , SOH^+ interference on Ti isotopes

HR-ICP-MS effective technique, but expensive




Determination of Ti in biological samples using ICP-MS

- Preliminary work started to measure titanium in hip samples, via serum samples
- Three modes compared:- He KED, SQ NH₃ and TQ NH₃
- Aim: To test if TQ mode gives low enough LOQ to enable determination of the normal Ti levels in patient samples
- Lowest LOQ only possible with Ti isotope at m/z 48 (abundance 73.8%), but serum high in Ca (⁴⁸Ca interference)
- Solution: Use ammonia as the reaction gas to isolate m/z 48 Ti from Ca



Reaction of Ti with NH₃: how it works

- **Q1** – set to transmit Ti, potential interferences on the product ion (e.g. ¹¹⁴Cd) and lower mass interference precursors (e.g. ³¹P, ¹⁶O) rejected.
- **Q2** – filled with NH₃. Ti collides and generates a range of adducts including ⁴⁸Ti(NH₃)₃NH⁺ at mass 114
- **Q3** – set to transmit mass 114, other masses rejected.



Comparison of different ICP-MS modes for Ti analysis

Sample matrix - 1:10 diluted serum plus 1ppm Cd, all data in $\mu\text{g/L}$

Sample i.d.	He KED mode, on mass at ^{48}Ti	Ti SQ NH_3 mode, at mass 114	Ti TQ NH_3 mode, at mass 114	Ti reported value, measured at ^{47}Ti using HR-ICP-MS
Serum L-1	167	1800	6.64	6.8
Serum L-1	262	1850	6.38	6.8

^{48}Ca interference plus residual PO^+ etc.

Contribution from ^{114}Cd

Only TQ NH_3 mode is capable of providing the correct Ti result

Arsenic measurement in the presence of cobalt

- Determination of elemental impurities in Vitamin B12
- Vitamin B12 contains Co (approx. 4% (w/w))
- Elements to be measured As, Cd, Pb and Hg – the so-called ‘Big Four’ in pharmaceutical analysis
- Digest sample in nitric acid
- Run all elements in SQ-KED mode and also As in TQ-O₂ mode (as ⁷⁵As¹⁶O)




Performance in SQ and TQ modes

Concentration Vitamin B12	Signal at $m/z=59$ (SQ-KED) [CPS]	Signal at $m/z=75$ (SQ-KED) [CPS]	BEC in SQ-KED mode [$\text{ng}\cdot\text{g}^{-1}$]	Signal at $m/z=75$ (TQ-O ₂)	BEC in TQ-O ₂ mode [$\text{ng}\cdot\text{g}^{-1}$]	Spike recovery in TQ-O ₂ mode [%]
BLK	73	2	0.0008	4	0.0007	N/A
0.0001 $\text{mg}\cdot\text{mL}^{-1}$	202,455	13	0.003	9	0.001	100.1
0.001 $\text{mg}\cdot\text{mL}^{-1}$	2,174,144	88	0.02	10	0.001	99.5
0.01 $\text{mg}\cdot\text{mL}^{-1}$	24,003,087	852	0.21	8	0.001	101.8
0.1 $\text{mg}\cdot\text{mL}^{-1}$	243,093,619	8744	2.47	18	0.002	106.4

- SQ-KED mode elevated BEC due to CoO contribution that cannot be suppressed with He KED.
- TQ-O₂ mode - measure AsO at m/z 91 free from CoO interference
- Accurate spike recovery (1 ng/g As) achieved with increasing concentrations of Vitamin B12 in TQ-O₂ mode

iCAP TQ on mass measurement example

- Measurement of Yb in a Gd matrix
- Same number of isotopes
- Similar abundances
- 16 mass units apart




	Symbol	Mass	Abundance
<input type="checkbox"/>	152Gd	151.9198	0.20
<input type="checkbox"/>	154Gd	153.9209	2.18
<input type="checkbox"/>	155Gd	154.9226	14.80
<input type="checkbox"/>	156Gd	155.9221	20.47
<input type="checkbox"/>	157Gd	156.9240	15.65
<input type="checkbox"/>	158Gd	157.9241	24.84
<input type="checkbox"/>	160Gd	159.9271	21.86


	Symbol	Mass	Abundance
<input type="checkbox"/>	168Yb	167.9339	0.13
<input type="checkbox"/>	170Yb	169.9348	3.05
<input type="checkbox"/>	171Yb	170.9363	14.30
<input type="checkbox"/>	172Yb	171.9364	21.90
<input type="checkbox"/>	173Yb	172.9382	16.12
<input type="checkbox"/>	174Yb	173.9389	31.80
<input type="checkbox"/>	176Yb	175.9426	12.70

Yb in a Gd matrix SQ

- Calibration 0 – 5 ppb Yb in 10 ppm Gd – no gas
- Calibration 0 – 5 ppb Yb in 10 ppm Gd – KED
- NH_3 reacts with many of the polyatomic ions that interfere with the REE however NH_3 also reacts quickly with some REE
- Pr, Eu, Dy, Ho, Er, Tm and Yb are less reactive with NH_3




^{172}Yb , no gas mode



^{172}Yb , KED mode

Yb measurement in 10pm Gd – TQ NH₃ mode




- Sensitivity – 7100 cps/ppb
- BEC – 0.05 ppb
- IDL – 0.0001ppb

- Yb measured on mass at m/z 172
- NH₃ flow – 0.9 ml/min

Isotope ratio example - Pb in the presence of Hg and REE

- Pb/Pb dating in geochronology
- Non radiogenic isotope ^{204}Pb used to correct for lead naturally occurring
- ^{204}Pb used as reference isotope for which others are compared
- ^{204}Pb has direct spectral overlap from ^{204}Hg that could be present
- Difficult to resolve these peaks even with HR-ICP-MS
- Normally use mathematical equations which could introduce errors



Pb isotope ratio results with Hg added – SQ mode




Sample i.d	$^{204}\text{Pb}/^{208}\text{Pb}$
Average ratio	0.02671
1ppb Pb	0.02571
1ppb Pb + 5ppb Hg	0.40942
1ppb Pb + 10ppb Hg	0.82649
1ppb Pb + 20ppb Hg	1.61867

- Measure isotope ratios in SQ mode
- Solutions with increasing Hg concentration
- Isotope ratio increases with increasing m/z 204 intensity

Isotope ratio with Hg and Yb added – SQ NH₃ mode

Sample i.d.	²⁰⁴ Pb/ ²⁰⁸ Pb
Theoretical	0.02671
1ppb Pb	0.02571
1ppb Pb + 5ppb Hg	0.02572
1ppb Pb + 1ppm Yb	0.07960

- SQ mode using NH₃ in the QCell
- Hg reacts, so Pb interference free at m/z 204
- However, Yb forms NH₃ cluster that SQ mode cannot resolve



Isotope ratio with Hg and Yb added – TQ NH₃ mode

Sample i.d.	²⁰⁴ Pb/ ²⁰⁸ Pb
Theoretical	0.02671
1ppb Pb	0.02546
1ppb Pb + 5ppb Hg	0.02567
1ppb Pb + 10ppb Hg	0.02542
1ppb Pb + 20ppb Hg	0.02563
1ppb Pb + 1ppm Yb	0.02566

- Measurements repeated in TQ NH₃ mode
- Again, Hg reacts with NH₃, so Pb free from Hg interference at m/z 204
- Yb rejected by Q1 so cannot form NH₃ cluster interference on m/z 204
- **Accurate ²⁰⁴Pb/²⁰⁸Pb ratios obtained in TQ mode**

Standard mode (i.e. no cell gas) with SQ operation

He KED single quadrupole mode with cell pressurised with He and KED applied

TQ NH₃ / H₂ / O₂ triple quadrupole mode with CRC pressurised with reaction gas Q1 set to analyte mass and Q3 set to either analyte mass (on mass analysis) or product ion (mass shift analysis)

- Flexibility and usability of both single and triple quadrupole modes
 - Full multielemental analysis with dedicated TQ interference removal for difficult analytes and simple He KED mode for everything else **in one analytical run**

Fully integrated autosampler and autodilution solutions



Elemental Scientific prepFAST



CETAC SDX_{HPLD}



Fully integrated speciation (IC and LC) and laser solutions



Redefining TQ-ICP-MS – ChromControl for speciation



IC / LC / GC-ICP-MS with fully integrated software control
ChromControl plug-in based on Chromeleon™ 7.2 CDS



Questions?

thermo scientific

www.thermofisher.com/iCAPTQ

www.thermofisher.com/TQ-ICP-MS

Thermo Scientific
iCAP TQ ICP-MS

Redefining triple quadrupole ICP-MS
with unique ease of use

