

The new NanoDrop™ micro-volume spectrometer for DNA analytics


Erich Heiden
Area Manager EEMEA
Core UC; μ UV; NMR



- Nucleic Acid Workflow
- The NanoDrop concept
- NanoDrop Lite
- NanoDrop 8000
- NanoDrop One and Acclaro™ software
- NanoDrop for protein and cell growth
- NanoDrop in the workflow





How do life scientists use NanoDrop to evaluate samples?



- Sample evaluation is a key step at multiple points during many life science workflows.
- Samples result from extraction and purification steps. Residual contaminants can be introduced during this process.
- Samples are used in complex and costly downstream experiments. The success of the experiment depends in part to the quality of the sample.

Nucleic Acid concentration from A_{260} (UV measurement)

- Nucleic Acids are made of five bases, A, T, G, C and U
- A nucleic acid spectrum is a “blending” of all nucleotide spectra in the sample
 - *Beer- Lambert equation:* $A_{260} = c * \epsilon * b \Rightarrow c = (A_{260} * \epsilon) / b$
- Everything that absorbs around 260 nm will contribute to A_{260} value: dsDNA, ssDNA, free oligos, dNTPs, RNA, proteins(absorb at 280nm), some contaminants
- A_{280} and Spectra deliver information about sample purity & quality



- Contaminants (e.g. Protein , Phenol) are often found in NA samples
 - Extraction kits (Trizol)/ Phenol)
 - From a “dirty” sample
- Protein and Phenol absorb at 280nm, and can affect the A_{260} measurement
- **Good:**
 - $A_{260}/A_{280} > 1.8$ signals minimal amounts of protein and residual Phenol
- **Better:**
 - Add A_{260}/A_{230} between 1.8 – 2.2 to ensure no residual salts (e.g. Guanidine HCl)
 - Familiarity with the “right” spectrum shape can reveal excess salts and proteins
- **Best** (more on this later!)

The Thermo Scientific NanoDrop pedestal technology



- ✓ Pipette 1 -2 μL of sample directly onto the measurement surface



- ✓ Close the arm and sample is drawn into a liquid column
- ✓ Optical fibers transmit light from a Xenon flash lamp through sample column; sample absorbance is measured



- ✓ After measurement, sample is simply wiped away
- ✓ No Carryover



By moving the pedestals, absorbance is measured at multiple path-lengths to achieve a very wide dynamic range (2-27,500ng/ μL for dsDNA). **No dilutions needed**

The NanoDrop Product Line Today - Summary

Absorbance



- DNA/Protein quantitation
- Single sample
- Most basic

- DNA/Protein quantitation
- Other absorbance studies
- Single sample
- Full Spectra



- DNA/Protein quantitation
- Other absorbance studies
- Up to 8 samples at a time
- Full Spectra

Fluorescence



- FLR label quantitation
- Other fluorescence studies
- Single sample


The NanoDrop Product Line

NanoDrop-Lite



NanoDrop Lite UV-Vis Spectrophotometer

- **Compact design, fits on any bench, fits small budgets**
- **Designed for routine measurements of NA and proteins**
- Small color screen with built-in methods
- Uses LEDs to measure absorbance at 260 and 280 nm
 - Purified Nucleic Acid A260, A260/A280 purity ratio
 - Purified Proteins A260
- Fixed wavelengths, no spectral output
- Dynamic range 4 – 1,500 ng/μL dsDNA
- Last 500 measurements saved locally, save data to USB device
- Optional docking station style printer, prints cryogenic tube labels

Thermo SCIENTIFIC	dsDNA (Factor: 50)		
	#1		
	A260 (10 mm): 3.578		
	A260/A280: 2.11		
	178.9 ng/μl		
	29 Nov 2011 10:29:41 Blank: 29 Nov 2011 10:28:47		
Home	Blank	Print	Measure



Why choose the ND-Lite?

For routine measurements of purified DNA, RNA and proteins (A_{260} & A_{280})

- Augment Full-spectrum NanoDrop products
- Value-priced “first NanoDrop” when moving to Micro-volume UV quant of DNA





The NanoDrop Product Line

NanoDrop-8000




NanoDrop 8000 UV-Vis Spectrophotometer



- **Designed for labs with higher throughput needs**
 - **Biopharma, biobanking, CROs and CMOs**
- Full-spectrum UV-Vis absorbance (220-750nm)
- Dynamic range: 2 - 3,700 ng/mL (dsDNA)
0.15 - 100mg/mL (BSA)
- Single Sample or 8-Sample modes
- 8 samples simultaneously in 20 seconds
- PC software for instrument control and data analysis



NanoDrop 8000 Software




- Designed for life science labs,
- Software offers preconfigured modules
 - Nucleic Acid, Proteins, Colorimetric Protein
 - Build your own methods, including Factor, Standards, and Oligos
- Measure labeled nucleic acids and proteins using MicroArray and Proteins & Labels methods



NanoDrop 8000: Sample Position Illuminator for 96 well plates




- This guide lights the respective samples to be assayed next.
- The Sample Position Illuminator correlates directly to the color-coded samples on the software screen.
- The Column Auto-advance feature identifies the next samples awaiting measurement.
- Can analyze 96 samples with full spectral data in <6min



Why choose the NanoDrop-8000?

For medium-throughput measurements of small-volume DNA, RNA and proteins in plate-based workflows

- Higher-throughput with only 1-2 μ l
- Full-spectrum analysis
 - Know your Purity ratios at 230 nm and 280 nm
 - Measure NA, Proteins and labeled samples
- PC software is intuitive and flexible
 - User Methods
 - Report generation, save, and print



NanoDrop-3300



NanoDrop 3300 Fluorospectrometer

- Fluorescence measurements from 1-2 uL of sample


Strengths

- Minimal sample required, no waste
- Flexibility –can use multiple fluorescence kits/vendors
- “Out of the box” applications




NanoDrop 3300 System Overview



- Fluorescence captured at 90° angle
- Excitation by one of three fixed LEDs
 - UV 365nm
 - blue 470nm
 - white 460nm-650nm
- CCD Detector records emission wavelengths between 400nm-750nm
- 1-2 μL sample
- **Communication, power via USB**
- Fluorescence signal output into a custom software interface



NanoDrop 3300 Excitation Source Overview



ND 3300 Applications –Protocols available



Nucleic Acid

- PicoGreen® assay dsDNA
- RiboGreen® assay RNA
- Hoechst 33258 dsDNA
- Sybr Green I
- Quant-It™ DNA HS
- Quant-It™ DNA BR

Proteins

- Quant-iT™ Protein assay
- FluoroProfile
- Fluoraldehyde OPA
- Fluorescamine


Other Samples

- 4-Methyl Umbelliferone
- Cy3/Alexa555
- Cy5/Alexa647
- Fluorescein-FITC-FAM
- Quinine Sulfate
- Dylights- 405,488,549,633,649, 680

Fluorescence Profiler



- Method generation tool
 - Determines the fluorescence maxima
 - Suggests best excitation
 - Direct method saving into the method list
 - Able to save multiple methods from one measurement



- Nucleic acid quantitation
- Protein quantitation
- Molecules with inherent fluorescence (GFP)
- Fluorescent Probes

- Forensics
- FRET
- Molecular Beacons
- Immunochemistry
- Microgenomics
- Histocompatibility
- 4-Methyl Umbelliferone
- Cy3/Alexa555
- Cy5/Alexa647
- Fluorescein-FITC-FAM
- Quinine Sulfate
- Evitag-UV

Thermo SCIENTIFIC NanoDrop Products

Home Overview Products **Science** Support Order News Contact Us Download Software

Science

How It Works

Application Notes

Fluorescence

UV/Vis

Protocols

Publications

Citations

Fluorescence

Measuring Enzyme Activity Using the Thermo Scientific NanoDrop 3300 Fluorospectrometer

A High-throughput Method for Measurement of Glomerular Filtration Rate in Conscious Mice (JoVE Protocol Video)

Measuring Mouse GFR by FITC-Inulin using the Thermo Scientific NanoDrop 3300 Fluorospectrometer

Microvolume Quantification of Proteins by UV-Vis Absorbance or Fluorescence (Poster)

Using the NanoDrop 3300 Fluorospectrometer for Cancer Research Applications

Using the NanoDrop 3300 Fluorospectrometer For Micro-Volume Fluorescence/FRET Applications

Using the NanoDrop 3300 Fluorospectrometer For Nanoparticle/Quantum Dot Applications

Molecular Beacon Probes: Micro-Volume Fluorescence Measurement of HPLC Isolated Probes

© Copyright 2013 Thermo Fisher Scientific Inc. All rights reserved. [Visit Thermo Fisher Corporate Website](#) | [Privacy Statement](#) | [Terms and Conditions](#)

The NanoDrop Product Line

NEW, unique NanoDrop One/ One^C



Introducing: NanoDrop One and NanoDrop One^c



- Will replace the current NanoDrop 2000/2000c
- Offer unique features that address user pain points
- Maintain clear leadership position of NanoDrop products demonstrating
 - Commitment to **innovation**
 - Commitment to **customer satisfaction** with compelling, relevant features

The NanoDrop One Value Proposition

- **The NanoDrop promise**

- Fast and easy answers from microvolumes
- Now even better

- **Acclaro™ Sample Intelligence**

- Contaminant identification using *data analysis* algorithms, supports downstream success
- On-demand embedded technical support answers context-sensitive questions
- Embedded camera with digital image analysis ensures measurement integrity

- **A high quality experience**

- New modern design with high resolution touchscreen
- Local control software with workflow-centric design
- Enhanced connectivity, extended data management options




NanoDrop One: The NanoDrop promise enhanced




Fast and easy to use

- NanoDrop One - even *faster* with Auto-Measure
- NanoDrop One - even *easier* with local control and touch screen interface



Microvolume measurements

- No need for slides or cuvettes, pipette directly on pedestal
- No waste of precious sample, measurements with 1-2 μL



No need for dilutions

- NanoDrop One - five auto-range pathlengths 0.03-1.0 mm, more accurate for concentrated samples
- NanoDrop One - extended dynamic range 2-27,500 ng/ μL dsDNA, *better accuracy*

Five auto-adjust pathlengths



Beer-Lambert equation

$$A = c * \epsilon * b$$

Shorter pathlengths (**b**) >> higher concentrations (**c**) can be measured

NanoDrop One offers five pathlengths, 0.03 -1.0 mm

- Higher dynamic range, no dilutions needed
- More accurate results for concentrated samples
- Can measure up to 400mg/mL IgG

NanoDrop 1000
2 pathlength

0.2 mm
max conc.=3,750
ng/ μ L dsDNA



NanoDrop 2000/c
4 pathlengths


0.05 mm
max conc.=15,000
ng/ μ L dsDNA



NanoDrop One/One^c
5 pathlengths

0.03 mm
max conc.=27,500
ng/ μ L dsDNA

How do life scientists use UV-Vis to evaluate samples?



- Sample evaluation is a key step at multiple points during many life science workflows.
- Samples result from extraction and purification steps. Residual contaminants can be introduced during this process. Small yields, final sample is precious.
- Samples are used in complex and costly downstream experiments. The success of the experiment depends in part to the quality of the sample.

Pain points in life-science sample evaluation

Instrument not designed for me ✓

- UV-Vis spectrophotometers are not made for “my kind” of samples
- Require training, cleaning cuvettes, and manual calculations
- Time consuming

Unknown concentration ✓

- How do I know I am within dynamic range?
- Dilutions introduce errors

Samples are precious ✓

- I have to use half of my sample to measure its concentration
- I wish I could evaluate my samples more than once but I don't have enough

Sample quality is critical ✓

- How do I know if my sample is good enough to use in my next experiment?
- Contaminated or inaccurately quantified samples lead to experimental failures

VoC indicates customers are concerned about sample quality

Bioinformatics survey

- 90% of customers look at spectra
- Customers look at purity ratios
- Customers concerned about contaminants in NA samples

-

Protein in DNA (RNA) samples	83%
Ethanol	72%
Phenol	57%
Guanidine	37%
DNA in protein	83%



The NanoDrop One Value Proposition

- **The NanoDrop promise**

- Fast and easy answers from microvolumes
- Now even better

- **Acclaro™ Sample Intelligence**

- Contaminant identification using *data analysis* algorithms, supports downstream success
- On-demand embedded technical support answers context-sensitive questions
- Embedded camera with digital image analysis ensures measurement integrity

- **A high quality experience**

- New modern design with high resolution touchscreen
- Local control software with workflow-centric design
- Enhanced connectivity, extended data management options





- **Identifies contaminants** in sample using powerful data analysis algorithms



- **Ensures measurement integrity** with embedded sensor and digital image analysis that monitors for broken columns and bubbles




- **Delivers on-demand** technical support with information alerts and guided troubleshooting



Acclaro = “to make clear”

Contaminant identification



Acclaro flags dsDNA sample #4


Acclaro identifies:

Nucleic Acid contaminants:

- Protein, phenol, guanidine HCl, guanidinium SCN

Protein contaminants:



- Nucleic acids and phenol




Protein in the sample contributes to the A260 and inflates dsDNA concentration. Protein contribution is subtracted from the original value to give the corrected concentration.

Ensure measurement Integrity

- Image of the sample column collected by embedded sensor is digitally analyzed
- Broken column:
 - Cancel measurement
- Bubbles in sample:
 - the air in bubbles diffracts light, the spectrum quality is degraded and the calculated concentration is inaccurate




User Alert: Bubbles Detected



What do I do now???

Learn more with Acclaro...



The screenshot shows the 'Acclaro Sample Intelligence' screen. It features a 'Bubbles Detected' notification with a 'Learn more' button. Below the notification, the text reads: "These measurement results may be invalid. Possible causes: Bubbles in sample or blank, which can cause poor liquid column formation and diffract UV light. Both issues will impact reproducibility. Particulates in sample or blank, which can diffract UV light and affect reproducibility. Possible solutions: To remove bubbles, briefly centrifuge sample and/or blank at low speed. To remove reflective particles such as magnetic beads, try extracting the sample again taking care to avoid particle contamination. For DNA samples, heat to 63 °C for 10 minutes, then centrifuge at low speed. For other sample types, centrifuge at high speed to condense particulates before taking an aliquot." A 'Done' button is in the top right corner.

The screenshot shows the 'Acclaro Sample Intelligence' screen with a 'Bubble Error Spectrum' section. It includes a text description: "When the sample or blank contains bubbles, the resulting spectrum may be irregular and the calculated analyte concentrations are usually invalid." Below this is a smaller graph showing a spectrum with a peak at 260 nm. A data table for this spectrum shows values for '#', 'mg/mL', 'A280', and 'A260/A280':


#	mg/mL	A280	A260/A280
1	0.3891	0.39	2.12

The graph is captioned "Spectrum of sample that contains large number of bubbles". A 'Learn More' button is at the top right of the section. The background shows a 'Measure' button set to 'ON'.


The screenshot shows the 'Acclaro Sample Intelligence' screen with two 'Learn More' options: "Spectrum of Sample with Bubbles" and "Effects of Bubbles in Sample". A 'Done' button is in the top right corner. The background shows a 'Measure' button set to 'ON'.



- Based on the most common questions received by ND tech support
 - Educates user on the science with short video animations
 - Presents common causes of particular problem, recommends solutions
 - Guides through troubleshooting with support tools
- Deliver information on sample quality and provide guided troubleshooting
- Topics covered are:
 - Purity ratios: when A_{260}/A_{280} or A_{260}/A_{230} ratios are outside acceptable limits for a pure DNA or RNA sample
 - Purity ratios: when A_{260}/A_{280} is outside acceptable limits for protein sample (A_{280})
 - Effects of dirty pedestals
 - Blanking issues
 - Baseline correction




User Alert: Purity ratios unreliable




When DNA&RNA sample concentration is <math><20\text{ng}/\mu\text{L}</math> the purity ratios are not reliable

The screenshot shows the 'Acclaro Sample Intelligence' alert dialog box. The title is 'Acclaro Sample Intelligence' with a 'Done' button. The main heading is 'Purity Ratios Unreliable'. Below this is a 'Learn more' button. The text in the dialog reads: 'When the DNA or RNA sample concentration is approximately 20 ng/μL or less, the purity ratios are unreliable because the sample spectrum is fairly flat and the calculated purity ratios are close to zero.' The background shows the same software interface as the first screenshot.

User Alert: Purity Ratio Out of Range



What do I do now???



A260/A280 Purity Ratio Below Acceptable Limit

[Learn more](#)


The A260/A280 ratio for this sample is below 1.65 due to strong absorbance at 280 nm. For pure dsDNA and pure RNA this ratio is typically between 1.8 and 2.2.

Possible causes:

- Incorrect blanking solution
- Contaminate absorbance at 280 nm or lower, possibly from:
 - Residual extraction reagent (e.g., phenol)
 - Presence of proteins

Possible solutions:

- Measure a new blank using the same buffer solution used to resuspend the sample
- Extract the sample again taking care to avoid contamination, or contact extraction kit manufacturer for additional cleanup steps to optimize the extraction



A260/A280 Below 1.65 - Learn More

[Common Contaminants](#)

[video What is a Purity Ratio?](#)

[Reading Resources](#)



Take-aways:

- ✓ *Acclaro provides users with information about sample quality so they can make informed decisions on the suitability of the sample for use in downstream applications*
- ✓ *Sample quality is crucial for success in downstream experiments*
 - Successful experiments lead to publications, patents, PhDs, even tenure
 - Successful experiments avoid lengthy troubleshooting, costly delays, and having to repeat expensive procedures

NanoDrop One customers perform successful experiments

The NanoDrop One Value Proposition

- **The NanoDrop promise**

- Fast and easy answers from microvolumes
- Now even better

- **Acclaro™ Sample Intelligence**

- Contaminant identification using *data analysis* algorithms, supports downstream success
- On-demand embedded technical support answers context-sensitive questions
- Embedded camera with digital image analysis ensures measurement integrity

- **A high quality experience**

- New modern design with high resolution touchscreen
- Local control software with workflow-centric design
- Enhanced connectivity, extended data management options



Modern design, looks like no other

- Distinctive design promises a high quality experience
- Local Control reduces footprint, streamlines workflows
- High resolution touchscreen interface
 - tilts for good visibility whether user is standing or sitting
 - slides to fit the space and accommodate R-and L-handed users
 - is compatible with lab gloves
- Glass disk is lit when instrument is ON and ready to use
 - Flashes during measurement
- New: power switch



Design Patent Pending!


Run with Local Control, PC compatible

- Local Control
 - Streamlined workflow guides users through measurement
 - Applications-based software
 - Print data using Dymo printer
- PC software
 - For data viewing, analysis, saving, and printing
 - Can explore all the Acclaro support tools
 - User can create new custom methods on PC then run them on LC
- Enhanced connectivity, extended data management options
 - Wi-Fi connects you to network,
 - Bluetooth for mouse and keyboard
 - 3 USB-A ports






New protein capabilities

83% of our customers quantify proteins with a spectrophotometer at least once per week



- Protein A205 (new)
 - for proteins without Trp, Tyr
- Protein A280
 - Protein Editor in Settings - create and save your “special protein” (new)
- Proteins and Labels
- Protein colorimetric assays




Beyond NA and Proteins Applications



Software anticipated at FCS


- Nucleic Acids Applications
- Proteins Applications
- OD600 Applications
 - Enter conversion factor to convert OD600 value to #cells/mL(new)
- Custom Methods
 - Formula method without standard curve
- PC Viewer Software
- Ethernet tether NDOne to computer
- Sensor diagnostics

- Very easy to upgrade Local Control software
 - Download from our website on USB and update LC






Future releases coming by end of 2016

- UV-Vis “Scan” method
- UV peak find, 3-pt net
- Acclaro Blank evaluation 
- Acclaro Baseline correction 
- Custom Methods – standard curve
- Custom Methods - advanced features
- Kinetics (5 stages, three wavelengths)
- WiFi connectivity - send data to one or more computers; send data to network
- Translations – Spanish, Portuguese, Japanese, Chinese, German, French, Polish, English, Korean
- IQ/OQ



NanoDrop One Help System





The Nano Drop range

And now

Questions??

Thank you!

